



Site Preparation: User Supplied Material Checklist for SP-G2 DNA Isolation and DLS-G2 DNA Labeling Protocols

DOCUMENT NUMBER:

CG-00050

DOCUMENT REVISION:

C

EFFECTIVE DATE:

06/12/2024

Table of Contents

Revision History	3
Introduction	4
Site Details and Points of Contact	4
User Supplied Material and Equipment for SP-G2 DNA Isolation Protocols	5
Common to all Bionano SP-G2 DNA Isolation Protocols	5
SP-G2 Frozen Human Blood and Frozen Human BMA Protocol	6
SP-G2 Fresh Human BMA Protocol	6
SP-G2 Fresh Human Blood Protocol	7
SP-G2 Frozen Cell Pellet Protocol	7
SP-G2 Fresh Cell Pellet Protocol	7
User Supplied Material and Equipment for DLS-G2 DNA Labeling Protocol	8
Note on Disposal of Hazardous Waste	9
Acknowledgement	10
Customer Representative name, job title (printed):	10
Bionano FAS name (printed):	10
Technical Assistance	11
Legal Notice	12
Patents	12
Trademarks	12

Revision History

REVISION	NOTES
A	Initial release.
B	Updated to include SP-G2 Fresh BMA and Europe Tech Support Contact info.
C	Removed references to Saphyr to make document generally applicable; minor corrections and formatting edits

Introduction

This completed checklist will confirm that a customer to be trained on the Bionano SP-G2 DNA Isolation workflows for Fresh and Frozen Human Blood, Fresh and Frozen Human Bone Marrow Aspirate (BMA), and Fresh and Frozen Cultured Cells, as well as the DLS-G2 DNA Labeling workflow, has all required and/or recommended user-supplied materials and equipment. The relevant lab personnel must go through the checklist, take inventory of their lab, and acquire any items not already present in their lab. Upon receipt of the items, the lab personnel will check the 'Present' box in each row corresponding to the item and return the completed checklist to the Bionano Field Application Scientist (FAS).

Please note that new user onsite training will be scheduled only upon the completion and submission of this signed form. Please list the wet lab contact(s), as well as any other contact(s) who will oversee this process, to facilitate communication.

Please complete only the sections relevant to the type(s) of starting material(s) users expect to be trained on. Should you have questions on any specific items in the list, or on which sections should be completed, please contact your Bionano FAS or email support@bionano.com to discuss.

Site Details and Points of Contact

Institution name	Address	e-mail address and phone number
------------------	---------	---------------------------------

Point of Contact Dept.	Name	E-mail address
------------------------	------	----------------

User Supplied Material and Equipment for SP-G2 DNA Isolation Protocols

Common to all Bionano SP-G2 DNA Isolation Protocols

Item	Supplier	Catalog	Item Present
-20 °C freezer	General lab supplier		<input type="checkbox"/>
4 °C refrigerator	General lab supplier		<input type="checkbox"/>
DynaMag™-2 Magnet*	Thermo Fisher	12321D	<input type="checkbox"/>
HulaMixer™ Sample Mixer*	Thermo Fisher	15920D	<input type="checkbox"/>
Microcentrifuge Tubes, 2 ml, Nuclease Free	Fisher Scientific	05-408-138	<input type="checkbox"/>
Microcentrifuge Tubes, 5 ml, Nuclease Free	Eppendorf	0030119487	<input type="checkbox"/>
Ethanol, 200 Proof, Molecular Biology Grade	Sigma-Aldrich	E7023	<input type="checkbox"/>
Isopropanol (IPA), ≥ 99.5%, Molecular Biology Grade	Fisher Scientific	AC327270010	<input type="checkbox"/>
Conical Centrifuge Tubes, 15 ml, PP	Fisher Scientific	05-539-12	<input type="checkbox"/>
Conical Centrifuge Tubes, 50 ml, PP	Fisher Scientific	14-432-22	<input type="checkbox"/>
Microcentrifuge with 1.5 ml Tube Rotor (Refrigeration at 4°C needed for Fresh Cells Protocol)	Eppendorf or Equivalent	5425/5425R	<input type="checkbox"/>
Minicentrifuge	Fisher Scientific or Equivalent	12-006-901	<input type="checkbox"/>
Ice Bucket and Ice	General Lab Supplier		<input type="checkbox"/>
Eppendorf ThermoMixer® C, with Eppendorf SmartBlock™ 1.5 mL (Alternatives: Heat block or Water bath set to 55°C)	Eppendorf or Equivalent	5382000023 and 5360000038	<input type="checkbox"/> <input type="checkbox"/>
Pointed Forceps	Electron Microscopy Sciences or Equivalent	78141-01	<input type="checkbox"/>
Wide-Bore Pipette Tips, Filtered, Aerosol, 200 µl	USA Scientific or Labcon ZAP or Equivalent	1011-8810 or 1152-965-008-9	<input type="checkbox"/>
Filtered Extra Long 1000 µl Tips, Sterile	VWR or Equivalent	76322-154	<input type="checkbox"/>
Pipettes (10, 20, 200, and 1,000 µl) and Nuclease Free, Filtered Pipette Tips	General Lab Supplier		<input type="checkbox"/>

Item	Supplier	Catalog	Item Present
Aluminum Cooling Block for 1.5 ml and 2 ml (optional)	Sigma Aldrich or Equivalent	Z743497	<input type="checkbox"/>
Bath Sonicator (optional)	General Lab Supplier		<input type="checkbox"/>
Benchtop Vortexer	VWR or Equivalent	10153-838	<input type="checkbox"/>
Positive-Displacement Pipette MR-10 (recommended)**	Rainin or Equivalent	17008575	<input type="checkbox"/>
Pipette Tips, 10 ul, C-10 for Positive Displacement Pipette (recommended) **	Rainin or Equivalent	17008604	<input type="checkbox"/>
Qubit™ Assay Tubes*	Thermo Fisher	Q32856	<input type="checkbox"/>
Qubit™ 4 Fluorometer*	Thermo Fisher	Q33238	<input type="checkbox"/>
Qubit™ BR (Broad Range) dsDNA Assay Kit*	Thermo Fisher	Q32853	<input type="checkbox"/>

*(Items in blue) We strongly encourage not substituting this equipment and getting the exact item recommended, as using this will give the most successful outcome from the workflow.

**A positive displacement pipette and tips work using a disposable piston that moves within a plastic capillary. This pipette can be highly effective in accurately pipetting and dispensing small volumes of viscous liquids such as the ultra-high molecular weight gDNA generated from Bionano's SP-G2 DNA Isolation Protocols.

SP-G2 Frozen Human Blood and Frozen Human BMA Protocol

Item	Supplier	Catalog	Item Present
HemoCue® WBC Analyzer*	Fisher Scientific (for US) or Distributor (outside US)	22-601-017	<input type="checkbox"/>
HemoCue® Microcuvettes*	Fisher Scientific	22-601-018	<input type="checkbox"/>
Water Bath, 37°C	General Lab Supplier		<input type="checkbox"/>
Parafilm	General Lab Supplier		<input type="checkbox"/>

*(Items in blue) We strongly encourage not substituting this equipment and getting the exact item recommended, as using this will give the most successful outcome from the workflow.

SP-G2 Fresh Human BMA Protocol

Item	Supplier	Catalog	Item Present
HemoCue® WBC Analyzer*	Fisher Scientific (for US) or Distributor (outside US)	22-601-017	<input type="checkbox"/>
HemoCue® Microcuvettes*	Fisher Scientific	22-601-018	<input type="checkbox"/>
Parafilm	General Lab Supplier		<input type="checkbox"/>

*(Items in blue) We strongly encourage not substituting this equipment and getting the exact item recommended, as using this will give the most successful outcome from the workflow.

SP-G2 Fresh Human Blood Protocol

Item	Supplier	Catalog	Item Present
HemoCue WBC Analyzer*	Fisher Scientific (for US) or Distributor (<u>outside US</u>)	22-601-017	<input type="checkbox"/>
HemoCue Microcuvettes*	Fisher Scientific	22-601-018	<input type="checkbox"/>
Test Tube Rocker	VWR or Equivalent	10159-754	<input type="checkbox"/>
Parafilm	General Lab Supplier		<input type="checkbox"/>

*(Items in blue) We strongly encourage not substituting this equipment and getting the exact item recommended, as using these will give the most successful outcome from the workflow.

SP-G2 Frozen Cell Pellet Protocol

Item	Supplier	Catalog	Item Present
Water Bath, 37°C	General Lab Supplier		<input type="checkbox"/>

SP-G2 Fresh Cell Pellet Protocol

Item	Supplier	Catalog	Item Present
Hemocytometer & Phase Contrast Microscope or Automated Cell Counter	General Lab Supplier		<input type="checkbox"/>
Refrigerated Centrifuge with Swinging Bucket Rotor for 15 ml Conical Tubes	Eppendorf or Equivalent	5804R	<input type="checkbox"/>
Cell Culture Media	General Lab Supplier		<input type="checkbox"/>
Serological pipettes	General Lab Supplier		<input type="checkbox"/>
Water Bath, 37°C	General Lab Supplier		<input type="checkbox"/>

User Supplied Material and Equipment for DLS-G2 DNA Labeling Protocol

Item	Supplier	Catalog	Item Present
HulaMixer Sample Mixer	Thermo Fisher	15920D	<input type="checkbox"/>
Thermocycler with heated lid*	General Lab Supplier		<input type="checkbox"/>
PCR tubes, 0.2 ml, thin-walled, flat cap, nuclease-free	Thermo Fisher or equivalent	AM12225	<input type="checkbox"/>
Microcentrifuge tubes, 0.5 ml, Amber, nuclease-free	USA Scientific or equivalent	1605-0007	<input type="checkbox"/>
Pipette tips, wide bore, filtered, 200 µl	USA Scientific or Labcon ZAP or Equiv.	1011-8810, 1152-965-008-9	<input type="checkbox"/>
Pipette tips, standard, filtered; 2, 10, 20 and 200 µl	General lab supplier		<input type="checkbox"/>
-20°C benchtop enzyme cooler	VWR or equivalent	414004-286	<input type="checkbox"/>
4°C aluminum cooling tube block	Sigma Aldrich or equivalent	Z740270	<input type="checkbox"/>
Forceps, pointed and curved	Electron Microscopy Sciences or Equiv.	78141-01	<input type="checkbox"/>
Pipettes (2, 10, 20 and 200 µl)	General Lab Supplier		<input type="checkbox"/>
Ice bucket and Ice	General Lab Supplier		<input type="checkbox"/>
Benchtop Vortexer	VWR or Equivalent	10153-838	<input type="checkbox"/>
Minicentrifuge	Fisher Scientific or Equivalent	12-006-901	<input type="checkbox"/>
Qubit 4 Fluorometer	Thermo Fisher	Q33238	<input type="checkbox"/>
Qubit® Assay Tubes	Thermo Fisher	Q32856	<input type="checkbox"/>
Qubit® HS (High Sensitivity) dsDNA Assay Kit	Thermo Fisher	Q32851	<input type="checkbox"/>
Bath sonicator (optional)	Branson or equivalent	CPX 952-119R	<input type="checkbox"/>
Positive-displacement pipette MR-10 (recommended)**	Rainin or equivalent	17008575	<input type="checkbox"/>
Pipette tips, 10 µl, C-10 for pos. displacement (recommended)**	Rainin or equivalent	17008604	<input type="checkbox"/>

*(Items in blue) We strongly encourage not substituting this equipment and getting the exact item recommended, as using these will give the most successful outcome from the workflow.

*Thermocyclers with either lid temperature regulation or the option to turn lid heating on/off are acceptable for our workflow.

**A positive displacement pipette and tips work using a disposable piston that moves within a plastic capillary. This pipette can be highly effective in accurately pipetting and dispensing small volumes of viscous liquids such as the ultra-high molecular weight gDNA generated from Bionano's SP-G2 DNA Isolation Protocols.

Note on Disposal of Hazardous Waste

Buffers Digestion Enhancer, LBB, and WB1 provided in Bionano's SP-G2 DNA Isolation kits contain guanidine hydrochloride (GuHCl). GuHCl is harmful if swallowed or inhaled and causes skin and eye irritation. DO NOT mix with bleach or acidic reagents. Liquid waste containing GuHCl should be safely decontaminated with a quaternary ammonium disinfectant before disposal in a hazardous waste stream. We recommend following local environmental, health and safety regulations for decontamination and disposal of all solutions mixed with GuHCl.

Bionano's DNA Isolation Protocols also generate biological waste. Please follow your institution's guidelines for disposal of such biological and/or infectious waste.

Acknowledgement

By signing below, I declare that the information provided is accurate, in fact and form, to the best of my knowledge.

I understand that the success of the onsite training on the Bionano workflow depends on the possession of all required user-supplied materials and equipment.

Finally, I acknowledge that I have read the note on disposal of hazardous waste.

Customer Representative name, job title (printed):

Signature:

Date:

Comments:

Bionano FAS name (printed):

Signature:

Date:

Comments:

Technical Assistance

For technical assistance, contact Bionano Technical Support.

You can retrieve documentation on Bionano products, SDSes, certificates of analysis, frequently asked questions, and other related documents from the Support website or by request through e-mail and telephone.

TYPE	CONTACT
Email	support@bionano.com
Phone	Hours of Operation: Monday through Friday, 9:00 a.m. to 5:00 p.m., PT US: +1 (858) 888-7663 Monday through Friday, 9:00 a.m. to 5:00 p.m., CET UK: +44 115 654 8660 France: +33 5 37 10 00 77 Belgium: +32 10 39 71 00
Website	www.bionano.com/support
Address	Bionano, Inc. 9540 Towne Centre Drive, Suite 100 San Diego, CA 92121

Legal Notice

For Research Use Only. Not for use in diagnostic procedures.

This material is protected by United States Copyright Law and International Treaties. Unauthorized use of this material is prohibited. No part of the publication may be copied, reproduced, distributed, translated, reverse-engineered or transmitted in any form or by any media, or by any means, whether now known or unknown, without the express prior permission in writing from Bionano Genomics, Inc. Copying, under the law, includes translating into another language or format. The technical data contained herein is intended for ultimate destinations permitted by U.S. law. Diversion contrary to U. S. law prohibited. This publication represents the latest information available at the time of release. Due to continuous efforts to improve the product, technical changes may occur that are not reflected in this document. Bionano Genomics, Inc. reserves the right to make changes in specifications and other information contained in this publication at any time and without prior notice. Please contact Bionano Genomics, Inc. Customer Support for the latest information.

BIONANO GENOMICS, INC. DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. TO THE FULLEST EXTENT ALLOWED BY LAW, IN NO EVENT SHALL BIONANO GENOMICS, INC. BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT BIONANO GENOMICS, INC. IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

Patents

Products of Bionano Genomics® may be covered by one or more U.S. or foreign patents.

Trademarks

The Bionano logo and names of Bionano products or services are registered trademarks or trademarks owned by Bionano Genomics, Inc. (“Bionano”) in the United States and certain other countries.

Bionano™, Bionano Genomics®, Saphyr®, Saphyr Chip®, Bionano Access™, Bionano Solve™, VIA™ software, Stratys™, Stratys™ Compute, Stratys™ Chip, and Bionano EnFocus™ are trademarks of Bionano Genomics, Inc. All other trademarks are the sole property of their respective owners.

No license to use any trademarks of Bionano is given or implied. Users are not permitted to use these trademarks without the prior written consent of Bionano. The use of these trademarks or any other materials, except as permitted herein, is expressly prohibited and may be in violation of federal or other applicable laws.

© Copyright 2024 Bionano Genomics, Inc. All rights reserved.