



Saphyr Site Preparation: User Supplied Material Checklist for Plant DNA Isolation Protocols

DOCUMENT NUMBER:

CG-00052

DOCUMENT REVISION:

A

Effective Date:

07/19/2023

Table of Contents

Legal Notice	3
Revision History	4
Introduction	5
User Supplied Material and Equipment for Plant DNA Isolation Protocols	6
Acknowledgement	9
Technical Assistance	10

Legal Notice

For Research Use Only. Not for use in diagnostic procedures.

This material is protected by United States Copyright Law and International Treaties. Unauthorized use of this material is prohibited. No part of the publication may be copied, reproduced, distributed, translated, reverse-engineered or transmitted in any form or by any media, or by any means, whether now known or unknown, without the express prior permission in writing from Bionano Genomics, Inc. Copying, under the law, includes translating into another language or format. The technical data contained herein is intended for ultimate destinations permitted by U.S. law. Diversion contrary to U. S. law prohibited. This publication represents the latest information available at the time of release. Due to continuous efforts to improve the product, technical changes may occur that are not reflected in this document. Bionano Genomics, Inc. reserves the right to make changes in specifications and other information contained in this publication at any time and without prior notice. Please contact Bionano Genomics, Inc. Customer Support for the latest information.

BIONANO GENOMICS, INC. DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. TO THE FULLEST EXTENT ALLOWED BY LAW, IN NO EVENT SHALL BIONANO GENOMICS, INC. BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT BIONANO GENOMICS, INC. IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

Patents

Products of Bionano Genomics® may be covered by one or more U.S. or foreign patents.

Trademarks

The Bionano logo and names of Bionano products or services are registered trademarks or trademarks owned by Bionano Genomics, Inc. ("Bionano") in the United States and certain other countries.

Bionano™, Bionano Genomics®, Saphyr®, Saphyr Chip®, Bionano Access™, and Bionano EnFocus™ are trademarks of Bionano Genomics, Inc. All other trademarks are the sole property of their respective owners.

No license to use any trademarks of Bionano is given or implied. Users are not permitted to use these trademarks without the prior written consent of Bionano. The use of these trademarks or any other materials, except as permitted herein, is expressly prohibited and may be in violation of federal or other applicable laws.

© Copyright 2023 Bionano Genomics, Inc. All rights reserved.

Revision History

REVISION	NOTES
A	Initial release.

Introduction

This checklist is intended to confirm that a customer looking to get trained on the Bionano Plant DNA Isolation workflows has all required and/or recommended user-supplied material and equipment. The wet lab contact at the site is required to go through the checklist and order/procure any outstanding items. Upon receipt of outstanding items, user is required to check the 'Present' box in each row corresponding to the relevant sample type/protocol and return the completed checklist to their Bionano Field Application Scientist.

Please note that new customer onsite training will be scheduled only upon the completion and submission of this signed form. Please list the wet lab contact(s), as well as any other contact(s) who will oversee this process to facilitate communication.

Please complete only the sections relevant to the type(s) of starting material(s) that you expect to be trained on. Should you have questions on a specific item(s) in the list, or on which sections you should complete, please contact your Bionano Field Application Scientist or support@bionano.com to discuss.

Site Details and Points of Contact

Institution name	Address	e-mail address & phone number

Point of Contact Dept.	Name	E-mail address

User Supplied Material and Equipment for Plant DNA Isolation Protocols

Common to all Bionano Plant DNA Isolation Protocols

Item	Supplier	Catalog	Item Present
Sample Handling and Plug Lysis – Day 1			
Formaldehyde 36.5%-38% in H ₂ O	Sigma-Aldrich or equivalent	F8775	<input type="checkbox"/>
Spermidine trihydrochloride	Sigma-Aldrich or equivalent	S2501, 85578 or 85580	<input type="checkbox"/>
Spermine tetrahydrochloride	Sigma-Aldrich or equivalent	S1141, 85605 or S2876	<input type="checkbox"/>
β-mercaptoethanol (βME)	Sigma-Aldrich or equivalent	M6250	<input type="checkbox"/>
40 μm and 100 μm cell strainers	VWR or equivalent	21008-949 and 21008-950	<input type="checkbox"/>
50 ml conical tube(s)	VWR or equivalent	21008-951	<input type="checkbox"/>
15 ml conical tube(s)	VWR or equivalent	21008-918	<input type="checkbox"/>
Petri dish (square, 12 x 12 cm)	Fisher or equivalent	12-565-224	<input type="checkbox"/>
Razor blades	General Lab Supplier		<input type="checkbox"/>
Plastic spatula	Corning catalog or equivalent	3005	<input type="checkbox"/>
Small soft-bristled paint brush	General Supplier		<input type="checkbox"/>
Heat blocks or water baths set to 43°C and 70°C*	General Lab Supplier		<input type="checkbox"/>
Thermomixer C	Eppendorf or equivalent	5382000023	<input type="checkbox"/>
50 ml SmartBlock	Eppendorf or equivalent	5365000028	<input type="checkbox"/>
Proteinase K enzyme, Puregene, 5 ml	QIAGEN	158146	<input type="checkbox"/>
RNase A solution, 5 ml	QIAGEN or equivalent	158156	<input type="checkbox"/>
Green screened caps for 50 ml conicals	Bio-Rad or equivalent	1703711	<input type="checkbox"/>
Agarose plug molds and plug mold plungers	Bio-Rad or equivalent	1703713	<input type="checkbox"/>
Agarose, 2% CleanCut Low Melting Point	Bio-Rad or equivalent	1703594	<input type="checkbox"/>
Filtered wide-bore pipette tip, 200 μl	VWR or Rainin equivalent	46620-642	<input type="checkbox"/>
DNA Recovery – Day 2			
Metal spatula	VWR or equivalent	82027-530	<input type="checkbox"/>
Orbital shaker, 12 in x 12 in Platform	Cole-Parmer or equivalent	EW-51820-30	<input type="checkbox"/>
Agarase (0.5 unit/μl)	Thermo Fisher	EO0461	<input type="checkbox"/>
0.1 μm Dialysis membrane, 47 mm diameter	Millipore or equivalent	VCWP04700	<input type="checkbox"/>
Petri dish, 60 x 15 mm	VWR or equivalent	28384-092	<input type="checkbox"/>
TE Buffer, pH 8.0	Thermo Fisher Scientific or equivalent	AM9849	<input type="checkbox"/>
Non-filtered, beveled pipette tip, 200 μl	USA Scientific or equivalent	1111-1810	<input type="checkbox"/>
DNA Quantification – Day 3			
Bath sonicator	VWR or equivalent	82027-530	<input type="checkbox"/>

Vortexer			
Qubit Fluorometer 3.0	Thermo Fisher	Q33216	
Qubit BR Assay Kit	Thermo Fisher	Q32850	
Qubit Assay Tubes	Thermo Fisher	Q32856	

Plant Base Protocol

Item	Supplier	Catalog	Item Present
Rotor-stator homogenizer (+ disposable probes)	QIAGEN or equivalent	9001271 (990890)	<input type="checkbox"/>
Refrigerated centrifuge + swinging bucket – 4500 x g	Eppendorf™ equipped with A-4-44 Rotor or equivalent	5804R	<input type="checkbox"/>
High-strength clear centrifuge tubes, 15-16 ml, closures and adapter sleeves (for density gradient centrifugation only)	Kimble® 45600-15 Glass Class B 15 ml Centrifuge Tubes + caps, Nalgene® 3117-0160 16 ml Polycarbonate Centrifuge Tubes with closure 18, VWR cat. # 21010-911 or equivalent; Kimble® 45550-15 Centrifuge Tube Rubber Adapter Sleeve		<input type="checkbox"/>
Ice bucket and ice	General Lab Supplier		

Plant High Polysaccharide Protocol

Item	Supplier	Catalog	Item Present
Rotor-stator homogenizer (+ disposable probes)	QIAGEN or equivalent	9001271 (990890)	<input type="checkbox"/>
Refrigerated Centrifuge	Eppendorf or equivalent	5804R	<input type="checkbox"/>

Plant High Polyphenol Protocol

Item	Supplier	Catalog	Item Present
Rotor-stator homogenizer (+ disposable probes)	QIAGEN or equivalent	9001271 (990890)	<input type="checkbox"/>
Refrigerated centrifuge + swinging bucket – 4500 x g	Eppendorf™ equipped with A-4-44 Rotor or equivalent	5804R	<input type="checkbox"/>
High-strength clear centrifuge tubes, 15-16 ml, closures and adapter sleeves (for density gradient centrifugation only)	Kimble® 45600-15 Glass Class B 15 ml Centrifuge Tubes + caps, Nalgene® 3117-0160 16 ml Polycarbonate Centrifuge Tubes with		<input type="checkbox"/>

	closure 18, VWR cat. # 21010-911 or equivalent; Kimble® 45550-15 Centrifuge Tube Rubber Adapter Sleeve		
--	--	--	--

Liquid Nitrogen Protocol

Item	Supplier	Catalog	Item Present
Refrigerated centrifuge + swinging bucket – 4500 x g	Eppendorf™ equipped with A-4-44 Rotor or equivalent	5804R	<input type="checkbox"/>
High-strength clear centrifuge tubes, 15-16 ml, closures and adapter sleeves (for density gradient centrifugation only)	Kimble® 45600-15 Glass Class B 15 ml Centrifuge Tubes + caps, Nalgene® 3117-0160 16 ml Polycarbonate Centrifuge Tubes with closure 18, VWR cat. # 21010-911 or equivalent; Kimble® 45550-15 Centrifuge Tube Rubber Adapter Sleeve		<input type="checkbox"/>
Liquid Nitrogen	General Lab Supplier		<input type="checkbox"/>
Mortar and Pestle	VWR or equivalent	89038-148 and 89038-164	<input type="checkbox"/>

Acknowledgement

By signing below, I declare that the information provided is accurate, in fact and form, to the best of my knowledge.

I understand that the success of the onsite training on the Bionano workflow depends on the possession of all required user-supplied material and equipment. .

Customer Representative name, job title (printed):

Signature:

Date:

Comments:

Bionano FAS name (printed):

Signature:

Date:

Comments:

Technical Assistance

For technical assistance, contact Bionano Genomics Technical Support.

You can retrieve documentation on Bionano products, SDS's, certificates of analysis, frequently asked questions, and other related documents from the Support website or by request through e-mail and telephone.

TYPE	CONTACT
Email	support@bionano.com
Phone	Hours of Operation: Monday through Friday, 9:00 a.m. to 5:00 p.m., PST US: +1 (858) 888-7663
Website	www.bionano.com/support
Address	Bionano Genomics, Inc. 9540 Towne Centre Drive, Suite 100 San Diego, CA 92121
