



# **Bionano VIA™ MultiScale BAM Reference Builder 4.0 Release Notes**

DOCUMENT NUMBER:

RNOTE-00009

DOCUMENT REVISION:

A

Effective Date:

06/27/2023

## Table of Contents

<b>Legal Notice</b>	<b>3</b>
Patents	3
Trademarks	3
<b>Revision History</b>	<b>4</b>
<b>New Features in Version 4.0</b>	<b>5</b>
<b>Bug Fixes in Version 4.0</b>	<b>6</b>
<b>Technical Assistance</b>	<b>7</b>

---

## Legal Notice

### **For Research Use Only. Not for use in diagnostic procedures.**

This material is protected by United States Copyright Law and International Treaties. Unauthorized use of this material is prohibited. No part of the publication may be copied, reproduced, distributed, translated, reverse-engineered or transmitted in any form or by any media, or by any means, whether now known or unknown, without the express prior permission in writing from Bionano Genomics, Inc. Copying, under the law, includes translating into another language or format. The technical data contained herein is intended for ultimate destinations permitted by U.S. law. Diversion contrary to U. S. law prohibited. This publication represents the latest information available at the time of release. Due to continuous efforts to improve the product, technical changes may occur that are not reflected in this document. Bionano Genomics, Inc. reserves the right to make changes in specifications and other information contained in this publication at any time and without prior notice. Please contact Bionano Genomics, Inc. Customer Support for the latest information.

BIONANO GENOMICS, INC. DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. TO THE FULLEST EXTENT ALLOWED BY LAW, IN NO EVENT SHALL BIONANO GENOMICS, INC. BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT BIONANO GENOMICS, INC. IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

### **Patents**

Products of Bionano Genomics® may be covered by one or more U.S. or foreign patents.

### **Trademarks**

The Bionano logo and names of Bionano products or services are registered trademarks or trademarks owned by Bionano Genomics, Inc. (“Bionano”) in the United States and certain other countries.

Bionano™, Bionano Genomics®, Saphyr®, Saphyr Chip®, Bionano Access™, and Bionano EnFocus™ are trademarks of Bionano Genomics, Inc. All other trademarks are the sole property of their respective owners.

No license to use any trademarks of Bionano is given or implied. Users are not permitted to use these trademarks without the prior written consent of Bionano. The use of these trademarks or any other materials, except as permitted herein, is expressly prohibited and may be in violation of federal or other applicable laws.

© Copyright 2023 Bionano Genomics, Inc. All rights reserved.

## Revision History

REVISION	NOTES
A	Initial release.

---

## New Features in Version 4.0

1. Support for samples containing ChrM.
2. Command line utility to convert the binary msr file to a tab delimited text file (.tsv or .tsv.gz). Command to run: `java -jar msbam-reference-export.jar`.
3. Calculation of read depth now includes clipped reads and improper pair reads.
4. A new option to either use or not use BAM QC measures to select files for reference creation.
  - a. A checkbox selection **Exclude Input Files based on QC** has been added to the UI
  - b. If used, samples that fail (above threshold of 20%) will not be used. QC metric used to exclude files:
    - i. BaseQ<30
    - ii. MapQ<30
    - iii. Invalid Reads
    - iv. Nonusable Reads
  - c. The reference file header stores the value of the parameter along with the metrics and cut-off thresholds

## Bug Fixes in Version 4.0

These are the bugs fixed in 4.0 that existed in previous version(s) and were not included in an incremental build release.

1. When re-building a reference file using some BAM files and some raw depth files, all intermediate depth files get deleted incorrectly and result in a “file not found error.” Now, only the raw depth file associated with the selected BAM file is deleted before re-building the reference.

## Technical Assistance

For technical assistance, contact Bionano Genomics Technical Support.

You can retrieve documentation on Bionano products, SDS's, certificates of analysis, frequently asked questions, and other related documents from the Support website or by request through e-mail and telephone.

TYPE	CONTACT
<b>Email</b>	<a href="mailto:support@bionano.com">support@bionano.com</a>
<b>Phone</b>	Hours of Operation: Monday through Friday, 9:00 a.m. to 5:00 p.m., PST US: +1 (858) 888-7663
<b>Website</b>	<a href="http://www.bionano.com/support">www.bionano.com/support</a>
<b>Address</b>	Bionano Genomics, Inc. 9540 Towne Centre Drive, Suite 100 San Diego, CA 92121