



Bionano Access Installation Guide

DOCUMENT NUMBER:

CG-30170

DOCUMENT REVISION:

AH

Effective Date:

07/31/2023

Table of Contents

Legal Notice	4
Patents	4
Trademarks	4
Revision History	5
Introduction	6
Compatibility	7
System Requirements	8
Access Updates	9
For Windows	9
For Mac	9
For Linux (CentOS/Redhat)	9
Backup	10
Rollback	10
Prerequisites	11
For Windows	11
For Macintosh	11
For CentOS 7 and RHEL 7	11
Windows Installation	14
Installation Steps	14
Linux Installation	15
Macintosh Installation	16
Configuration	17

Example access.txt Configuration File	17
Configuration Settings	17
Starting the Web Server	19
User Accounts	21
Default Account	21
Technical Assistance	22

Legal Notice

For Research Use Only. Not for use in diagnostic procedures.

This material is protected by United States Copyright Law and International Treaties. Unauthorized use of this material is prohibited. No part of the publication may be copied, reproduced, distributed, translated, reverse-engineered or transmitted in any form or by any media, or by any means, whether now known or unknown, without the express prior permission in writing from Bionano Genomics, Inc. Copying, under the law, includes translating into another language or format. The technical data contained herein is intended for ultimate destinations permitted by U.S. law. Diversion contrary to U. S. law prohibited. This publication represents the latest information available at the time of release. Due to continuous efforts to improve the product, technical changes may occur that are not reflected in this document. Bionano Genomics, Inc. reserves the right to make changes in specifications and other information contained in this publication at any time and without prior notice. Please contact Bionano Genomics, Inc. Customer Support for the latest information.

BIONANO GENOMICS, INC. DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. TO THE FULLEST EXTENT ALLOWED BY LAW, IN NO EVENT SHALL BIONANO GENOMICS, INC. BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT BIONANO GENOMICS, INC. IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

Patents

Products of Bionano Genomics® may be covered by one or more U.S. or foreign patents.

Trademarks

The Bionano logo and names of Bionano products or services are registered trademarks or trademarks owned by Bionano Genomics, Inc. (“Bionano”) in the United States and certain other countries.

Bionano™, Bionano Genomics®, Saphyr®, Saphyr Chip®, Bionano Access™, VIA™ software, and Bionano EnFocus™ are trademarks of Bionano Genomics, Inc. All other trademarks are the sole property of their respective owners.

No license to use any trademarks of Bionano is given or implied. Users are not permitted to use these trademarks without the prior written consent of Bionano. The use of these trademarks or any other materials, except as permitted herein, is expressly prohibited and may be in violation of federal or other applicable laws.

© Copyright 2023 Bionano Genomics, Inc. All rights reserved.

Revision History

REVISION	NOTES
AB	Installation and dependencies updated for Access release 1.7.
AC	Updated hyperlinks for MAC computer installations
AD	Added instruction in prerequisite section.
AE	Adding step to install instruction to remove previous Access-Install folder.
AF	Added additional instructions in the prerequisite section.
AG	Updated Software download link
AH	Updates for Bionano Access 1.8 release and additional data backup guidance

Introduction

The purpose of this document is to detail how to install and configure Bionano Access®. This document will cover a variety of approaches to installation. It will cover multiple configurations including shared web server and standalone installations. It will also cover installation procedures for various operating systems including Linux, Mac, and Windows. The first part of this document will cover the process for a first-time installation.

On startup Bionano Access will check the configuration for Compute servers. During startup Access will connect and verify the proper Solve version is installed on any Compute server configured. If it cannot connect or verify the proper version of Bionano Solve, Bionano Access will not start. So, you will want to update the Bionano Solve version on your Compute servers before upgrading Bionano Access. The Compute servers will signal the Bionano Access server when job is complete. So be sure to assign the Bionano Access Server a static IP so the address the Compute servers use to communicate with the Access Server does not change.

A single shared instance of Bionano Access per organization is the recommended installation. It is not recommended to install Bionano Access on each user workstation or to have one instance for each instrument. This would defeat the collaborative design of Bionano Access. It is possible to run Bionano Access on a workstation in standalone mode if a Saphyr instrument or Compute server is available to visualize Bionano data.

Compatibility

Bionano Access is compatible with data generated by a Saphyr® instrument and bioinformatics output from the Bionano Solve pipeline. Bionano Access can be run on a virtual server if needed. Bionano Access is currently supported on the following operating systems:

- CentOS 7.x
- RHEL 7.x
- Macintosh (Catalina)
- Windows 10

System Requirements

16 GB of memory for standalone installations and 32 GB for shared web servers is recommended. For disk, estimate approximately 10 GB per assembly you plan to retain on the web server. The web server is very lightweight and will run easily on most systems. If you are noticing performance issues, contact support@bionano.com to help identify a configuration that will work best for the situation.

Access Updates

The prerequisite software packages (see Prerequisites section below) only need to be installed once. To update Access, perform the Access installation again (see section Windows Installation or Linux/Macintosh Installation). The Bionano Access installation can be run as many times as desired. It will not remove or overwrite any of your data, files, or configurations. You can determine the current version of Bionano Access installed by looking at the **Access Help** page (<http://<servername>:3005/Help/HelpIndex.html>).

For Windows

- Stop the Bionano Access web server.
- Go to the control panel, select programs and features.
- Uninstall Node.js.
- Run the installation for Node.js 12.22.10.1.

For Mac

- Stop the Bionano Access web server.
- Run the installation for Node.js 12.22.10
 - Apple may prevent the installation. Go to **System Preferences > Security & Privacy** to run installer.
- Upgrade npm
 - `sudo npm install npm@8.10.0 -g`
- install Java SE Development Kit 11.0.16.1 (May require user to create an Oracle account to download installer.)

For Linux (CentOS/Redhat)

- Stop the Bionano Access web server.
- Issue the following commands.
NOTE: User terminal may need to be restarted after issuing the curl command.

```
sudo yum -y update
sudo yum remove nodejs
sudo yum install https://rpm.nodesource.com/pub_12.x/el7/x86_64/nodejs-12.20.1-1nodesource.x86_64.rpm
sudo npm install npm@8.10.0 -g
sudo yum install java-11-openjdk-devel
```

Backup

Data backup is the responsibility of the end user facility and is strongly recommended on a recurring basis. Prior to updating the software, it is strongly recommended to conduct a backup of the data on the Bionano Access Server, including the databases for Access and VIA software. Administrator access for the server is required for database directory permissions, if not already established please reach out to Support@Bionano.com for assistance. Bionano Support may not set up database backup for the end user however guidance can be provided. For additional details, please refer to the *Saphyr Site Preparation Guide (CG-30145)*.

Rollback

If a backup system is in place, you can restore the backup to return the system to the previous version if needed. Bionano support should be contacted for assistance in the instance of needing to revert the BAS from a backup. Users should back up the environment on a regular basis.

Prerequisites

If the Bionano Access system will be integrated with a Compute cluster, Bionano Solve should be installed/updated in that environment first. Access will check that the correct version of the pipeline is installed before it starts.

Bionano Access requires PostgreSQL 12 (PostgreSQL 11 for Mac), Perl 5, Java 11, and Node.js 12. Links have been provided below to the appropriate installation instructions based on the operating system. Do not install different versions than what has been indicated. Be sure to perform these installations with administrator rights. During installation, Stack Builder is not required and use default paths and ports during install.

For Windows

PostgreSQL 12: <https://www.enterprisedb.com/downloads/postgres-postgresql-downloads>. Set password to “1rysvieW”, use port 5432, and remember to update the path.

Perl 5: <https://strawberryperl.com/>

Java 11: <https://www.java.com/en/download/>

Node.js 12: <https://nodejs.org/download/release/v12.22.10/node-v12.22.10-x64.msi>. Upgrade npm to 8.10.0 after node installation by executing ‘`npm install npm@8.10.0 -g`’ in command prompt.

For Macintosh

Python: <https://www.python.org/downloads/release/python-2715/>. Required by PostgreSQL for Mac.

PostgreSQL 11: https://sbp.enterprisedb.com/getfile.jsp?fileid=12569&_ga=2.216068630.1813352955.1591738038-1474724881.1591738038. Set password to “1rysvieW”, use port 5432, and do not install Stack Builder.

Perl 5: <https://www.perl.org/get.html>

Java 11: <https://www.oracle.com/java/technologies/javase/jdk11-archive-downloads.html#license-lightbox>. May require user to create an oracle account to download installer.

Node.js 12: <https://nodejs.org/download/release/v12.22.10/node-v12.22.10.pkg>. Upgrade npm to 8.10.0 after node installation by executing ‘`sudo npm install npm@8.10.0 -g`’ in terminal.

For CentOS 7 and RHEL 7

```
# update yum
sudo yum -y update

# install perl
sudo yum install -y perl
```

```
# install java
sudo yum install java-11-openjdk-devel

# install python3
sudo yum install -y python3

# install nodejs
curl --silent --location https://rpm.nodesource.com/setup_12.x | sudo bash -
sudo yum install -y nodejs

# update npm
sudo npm install npm@8.10.0 -g

# install postgresql
sudo yum install -y https://download.postgresql.org/pub/repos/yum/reporepms/EL-7-x86_64/pgdg-redhat-repo-
latest.noarch.rpm

sudo yum install -y postgresql12 postgresql-12-server postgresql12-contrib
sudo /usr/pgsqr-12/bin/postgresql-12-setup initdb
sudo cp /var/lib/pgsqr/12/data/pg_hba.conf /var/lib/pgsqr/12/data/pg_hba.conf.orig
sudo bash -c 'echo "local all all peer">/var/lib/pgsqr/12/data/pg_hba.conf'
sudo bash -c 'echo "host all all 127.0.0.1/32 md5" >>/var/lib/pgsqr/12/data/pg_hba.conf'
sudo bash -c 'echo "host all all ::1/128 md5" >>/var/lib/pgsqr/12/data/pg_hba.conf'

# start postgresql service
sudo systemctl start postgresql-12
sudo systemctl enable postgresql-12
sudo -i -u postgres psql -U postgres -d postgres -c "alter user postgres with password 'lrysv1ew!';"
sudo systemctl restart postgresql-12
```

RSA keys and keyless login to any servers that will be used for detection or assembly servers before running the installation the first time should be established. To generate new RSA keys, employ this command as the user:

```
ssh-keygen -t rsa
```

To set up keyless login for the Bionano account use this command where <<servername>> is the DNS name for that server or the IP address. Test it with an SSH connection afterwards. Users will not be prompted for a password if it worked.

```
ssh-copy-id <<servername>>
```

WARNING: Please make sure both PostgreSQL and Java are in the user path before attempting the Bionano Access Install. This can be tested by trying the 'psql' or 'java' commands at a command prompt. If they are not recognized the path will need to be updated. Example commands to update the path environment variable have been provided below, but they may vary depending on the exact system. **NOTE:** The Mac Catalina OS requires PATH to be added while in bash.

Mac:

```
PATH="/Library/PostgreSQL/11/bin:$PATH"  
export INSTALL4J_JAVA_HOME=/Library/Java/JavaVirtualMachines/jdk-11.0.16.1.jdk
```

Linux:

```
export PATH=$PATH:/usr/pgsql-12/bin/
```

Windows:

- Right click your computer in Windows Explorer and click **Properties**.
- Select **Advanced System Settings**.
- Click **Environment Variables** button.
- Under system variables (not user variables) select **Path**.
- Click **Edit** button.
- Add (if not already there) `C:\Program Files\PostgreSQL\12\bin`
- Add (if not already there) `C:\Program Files (x86)\Java\jdk1.8.0_161\bin`
- Click **Ok** (x3)

Windows Installation

64-bit Windows 10 or greater is supported. Follow the instructions below to complete a Windows installation.

Installation Steps

1. Make sure the system has Internet access.
2. Make sure prerequisites are installed and PostgreSQL and Java are added to PATH. (Follow page 14)
3. Make sure the active account has administrator rights.
4. Download the Bionano Access Windows Setup (`Access_Setup.exe`) from the Bionano Software Download page (<https://bionano.com/software-downloads/>)
5. Right click on the `Access_Setup.exe` file and select **Run as Administrator**.

Linux Installation

Most testing is performed using CentOS. As such CentOS is the recommended operating system. However, it is known that other organizations often prefer other flavors of Linux. Other versions of Linux are likely to be compatible as well but have not been validated for use at this time. The following instructions will guide how to install Bionano Access on a Linux system.

If not already applying security updates to your system on a regular basis we recommend running 'yum update' on your Linux operating system prior to installing the latest version of Bionano Access. This operation requires Internet access and will update the libraries on the system to the latest versions including various security updates. Be sure to reboot the system after applying the latest patches when your system is idle.

1. Verify the installation system has internet access.
2. Make sure prerequisites are installed and PostgreSQL is added to PATH. (Follow page 14)
3. Create a Bionano user account that has sudo privileges.
4. Download the Bionano Access Linux/Mac Installation (`Access-Install.tar.gz`) from the Bionano Software Downloads page (<https://bionano.com/software-downloads/>). Be sure to delete previous installation tar files before downloading the latest installation tar file.

5. Delete previous install folder:

```
rm -rf Access-Install
```

6. Extract the file to /home/bionano:

```
tar -zxvf Access-Install.tar.gz -C /home/bionano
```

7. Run installation script:

```
sudo /home/bionano/Access-Install/install.sh
```

Macintosh Installation

Mac installations are intended for standalone use only. We do not intend for Mac systems to function as shared web servers. You may be prompted for your Mac admin password during the installation. The install script will install nodejs if it is not already on your system.

1. Verify the installation system has internet access.
2. Make sure prerequisites are installed and PostgreSQL and Java are added to PATH. (Follow page 14).
3. Download the Bionano Access Linux/Mac Installation (`Access-Install.tar.gz`) from the Bionano Software Downloads page (<https://bionano.com/software-downloads/>) to your home directory.
4. Open a terminal window and go to the home directory (`cd ~/`).
5. Run this command: `bash`
6. Extract the install file in the home directory (`tar -xvf Access-Install.tar.gz`)
7. Go into the Access-Install directory (`cd ./Access-Install`)
8. Run this command: `./install_mac.sh`
9. Enter user password and PostgreSQL password ("`lrysviaw`") accordingly.

Configuration

A configuration file is used to specify the Bionano Access setup. This file tells Bionano Access where to find files and what Compute servers with which it should communicate. Specify the configuration file by name when starting the web server. If planning to run Bionano Access in standalone mode with no instruments or Compute servers it is not necessary to perform any configuration steps and can skip the balance of the Configuration section.

Bionano Access can be configured to work with a combination of Saphyr instruments and or Compute servers. Create an `access.txt` file in the `~/bionano/access/web/Server/Config` folder for the specific configuration. In the Config folder we provided various template files to use as a starting point. An `access.txt` file is not provided in the installation so that we do not overwrite your configuration during future updates.

Example `access.txt` Configuration File

```
{
    "psfDetectServer": "193.168.50.31",
    "assemblyServers": ["193.168.50.31"],
    "icsShare": "\\193.168.50.31\share",
    "webServerShare": "\\193.168.50.31\share",
    "irysSolveShare": "/home/bionano/access/share",
    "icsShareRelPath": "chips",
    "webServerShareJobs": "jobs",
    "webServerRelLocalFiles": "",
    "webServerPort": "3005",
    "solveServerUser": "bionano",
    "scriptDir": "/home/bionano/tools/access/1.0",
    "refAlignerDir": "tools/pipeline/1.0/REFALIGNER/1.0",
    "irysSolveListenIP": "10.10.1.13"
}
```

Configuration Settings

Adjust the following properties in the Access configuration file to establish proper connections with the Saphyr Instrument Control Software and user Compute servers.

Table 1. Connection properties

Property	Value	Example
psfDetectServer	This should be the IP address on the Saphyr Compute Server accessible to the Bionano Access Web Server. If you have more than one instrument it would be the primary Saphyr Compute server in your cluster. This value must be a static IP address.	"193.168.50.31"
irysSolveListenIP	This should be set to the IP address for the Bionano Access Web server that is accessible to all the Compute servers. In some cases, the web server may be dual homed so that the Compute servers can be within an isolated VLAN. In this case the web server may have more than one IP address. The Compute servers will curl http requests to the web server to update job status using this IP address. This value must be a static IP address.	"10.10.1.13"
icsShare	This is the path where the Saphyr Instrument Control Software will write images captured from the instrument for image detection. This path should be relative to the instrument control PC.	"\\\\193.168.50.31\\share"
webServerShare	This is the path to the share for images from the instrument and output for the pipeline. If the web server is running on a Windows machine this value should be the same as the icsShare property. If the web server is running on a Linux system it should be the same as the irysSolveShare property.	Windows: "\\\\193.168.50.31\\share" Linux: "/home/bionano/access/share"
assemblyServers	Input the static IP address of Compute node that can take SGE job submissions.	["193.168.50.31"]
webServerPort	The TCP port used by the Access application	3005
proxyServer	url of proxy server for internet access	http://192.168.48.21:3128
referenceListDisplay	Access ships with a variety of references for common studied species. It is possible to use this setting to limit references to human only. You can also limit references to only those that were provided with Bionano Access (not user uploaded).	"human;system"

Starting the Web Server

First open a command prompt and navigate to the following directory:

```
(Windows) c:\bionano\access\web\Server
(Linux) /home/bionano/access/web/Server
(MacOS) ~/access/web/Server
```

To start Bionano Access Web Server in standalone mode issue the following command:

```
node --max-old-space-size=16384 server default
```

To start Bionano Access Web Server with your access.txt configuration issue this command:

```
node --max-old-space-size=32768 server access
```

“node” is the command to start nodejs. The “--max-old-space-size” argument increases the V8 JavaScript engine memory use to be optimized for our larger data files. ‘server’ is our starting nodejs file and the final parameter (access or default) is the configuration file to use.

Add ‘https’ to the end of either command above if you install an ssl certificate on the server. It will need a key.pem and cert.pem in your access/web/Server directory. Please follow *How to enable HTTPS in Bionano Access (CG-30377)* to set up https for Access. We recommend using an SSL certificate for https communication to protect your data.

The web server should then be accessible from a browser using the IP address or DNS name for the webserver plus the port the server is running on. 3005 is our default port. Change the port for the web server in the configuration file by changing the ‘webServerPort’ setting. When enabling https the port will be the value of the webServerPort setting plus one. If using port 80 it is not necessary to add a port number to the URL, but there are known issues with curl responses from the Compute servers routing to port 80 on some networks due to security restrictions. If the Web Server’s IP address or DNS is unknown, contact the IT department for details, and if using standalone mode on a local machine such as a desktop or laptop, use ‘localhost’ in place of an IP address. Other users would require the IP address or DNS name instead of localhost to access the machine.

EX. <http://192.168.48.57:3005>, <http://bionano-227:3005/>

On CentOS and RedHat if closing your SSH session, the server will stop. Use the following commands to detach from the process and allow it to continue running after exit. When running this command the console output will be logged to the nohup.out file in the `.../bionano/access/web/Server` directory.

```
nohup node --max-old-space-size=32768 server access & disown
```

The first time Bionano Access 1.8 is started after upgrading from Access 1.6 or a previous version, several BED files will be jsonized. This process may take several minutes to complete. This conversion will only happen once and the Bionano Access webserver may not be responsive during this conversion.

NOTE: To start the web server in Linux, users can only use the “bionano” account which must have sudo privilege. Please do not start the access web server using the root account in Linux.

FIREWALL

There are a few security considerations with the Bionano Access web server. First, if firewall security is present, allow the following traffic from the Bionano Access web server. Generally, if running in standalone these considerations would not apply. The Compute On Demand endpoints will depend on the Compute region selected when activating the service. Contact support if you require the exact endpoints for these systems.

Table 2. Security System Information

Source System	Target System	Protocol
Bionano Access Web Server	psfDetectServer assemblyServers	SSH port 22
All Compute servers Saphyr Instrument Control PC Workstation Browsers	Bionano Access Web Server	HTTP/HTTPS (port is set as webServerPort in config file)
Compute On Demand	Regional Bionano Stratus Server Rescale Platform AWS	HTTP/HTTPS

User Accounts

When first accessing the web server, log in following the prompt. A default user account is provided to access the system. Once logged in, user accounts for all team members can be created.

WARNING: We suggest user the user create a new user account for Access usage. Do not use the default account for your analysis. We strongly recommend to disable or modify the default account once finished creating user accounts for the organization. Do not use shared accounts between users, which defeats how the notification scheme is designed to operate and it violates security guidelines for clinical environments. If changing the default account please notify Bionano Support or provide them with an account to use when servicing the system.

Default Account

Username: `bionanoadmin`

Password: `admin`

Technical Assistance

For technical assistance, contact Bionano Technical Support.

You can retrieve documentation on Bionano products, SDS's, certificates of analysis, frequently asked questions, and other related documents from the Support website or by request through e-mail and telephone.

TYPE	CONTACT
Email	support@bionano.com
Phone	Hours of Operation: Monday through Friday, 9:00 a.m. to 5:00 p.m., PST US: +1 (858) 888-7663
Website	www.bionano.com/support
Address	Bionano Genomics, Inc. 9540 Towne Centre Drive, Suite 100 San Diego, CA 92121