



Bionano Access™ Dashboard and Chip Metrics Guidelines

DOCUMENT NUMBER:

CG-30304

DOCUMENT REVISION:

F

Effective Date:

Aug-27-2026

Table of Contents

Revision History	3
Dashboard	4
Run Information	4
Analysis Graphs	5
Run Metrics Table	5
Run Troubleshooting	8
Chip Metrics Graphs	9
Technical Assistance	11
Legal Notice	12
Patents	12
Trademarks	12

Revision History

REVISION	NOTES
A	Initial release.
B	Document updated for Solve 3.7 release: - Update dashboard format and units - Add chip metrics graphs section
C	Updated document template. Updated Analysis Graphs and Run Troubleshooting sections.
D	Updated to include agnostic information intended for both Saphyr and Stratys instruments.
E	Minor updates including changing “QC Metrics” to “Run Metrics”
F	Hotfix update

Dashboard

This document provides guidelines on interpreting the quality metrics shown on the Bionano Access™ Dashboard. The Bionano Access Dashboard displays real-time graphic representations of ongoing Saphyr® or Stratys™ runs, updating throughput and quality metrics over time as the system scans. Dashboards from completed runs are archived within Bionano Access, so the user may refer back to data collection of a particular sample.

The dashboard consists of three main sections: run information, analysis graphs, and a run metrics table. It is a useful tool for evaluating the quality of the run—particularly how it proceeds over time, and how well molecule quality is sustained over long periods of data collection.

Run metrics generated in the dashboard are similar to those in the Molecule Quality Report (MQR). In **Figure 1** below, they are calculated as a weighted average per scan, updating over the course of the run. Once the run is complete, the dataset is imported into Bionano Access and a summary Molecule Quality Report is calculated. For the map rate and all related metrics (effective coverage, PLV, NLV), the dashboard uses molecules >150 kbp, whereas the MQR uses molecules > 150 kbp and minsites = 9. Label density is calculated using the entire population of molecules, not just molecules > 150 kbp.




Figure 1. Run metrics.

Run Information

- **Chip:** Chip serial number and lot number of a given run.
- **Run ID:** System-generated run identifier specific to a given run at a particular instrument.
- **Chip Name:** Chip name as entered by the user at Bionano Access.

- **Instrument:** Instrument serial number.
- **Min Length:** The cutoff of DNA length. Only DNA molecules with equal or longer length of the cutoff are counted in the Run Graphs.
- **Min Labels:** The cutoff of number of labels per molecule. Only DNA molecules with minimum or more number of labels are counted in the Run Graphs.
- **Start Time:** The time when the run was initiated. **NOTE:** this is relative to when data processing occurs in Access, not to when the chip is placed in the instrument.
- **End Time:** The time when the run was completed. **NOTE:** this is relative to when data processing completes in Access, not to when the chip is removed from the instrument.

Analysis Graphs

In **Figure 2**, the top panel displays DNA per scan (Gbp). This shows the total amount of DNA (>150 kbp) per scan in the flowcell as the run proceeds. The X and Y axes will scale to the data. The bottom panel shows the Map Rate (%), which is the percentage of molecules that map to the reference genome in each scan as the run proceeds. The X axis will scale as scans accumulate.

Users may click on either graph to zoom in, and click again to zoom out. If desired, a user can hide a flowcell line in the graphs by deselecting the checkbox in **Figure 3**.




Figure 2. Analysis graphs.

Run Metrics Table

The Run Metrics table is shown in **Figure 3** below. The average values are weighted averages based on the amount of DNA detected per scan.

Plot	Flowcell	Prep	Result	For Run ID: 2b178fb7-1e51-4a46-aae4-8cd63cbbfefe							All runs for chip: IV36UV6NPOOSFMGU	
				Avg N50 ≥150kbp	Avg N50 ≥20kbp	Avg Label Density per 100kbp	Avg Map Rate	Est Effective Coverage ≥150kbp	Avg PLV	Avg NLV	DNA Collected ≥150kbp	Cumulative DNA ≥150kbp
	1	Project: Chip QC Sample: 260701-0006153 Channel: Green 01 Enzyme: DLE-1 Reference: hg38_DLE1_0kb_labels_masked_YPARs.cmap		303 kbp	232 kbp	15.99	86.8%	177.362 X	8.1%	9.1%	Collected: 631.018 Gbp Target: 600.000 Gbp Percent: 105.170 %	631.018 Gbp
	2	Project: Chip QC Sample: 260701-0006153 Channel: Green 01 Enzyme: DLE-1 Reference: hg38_DLE1_0kb_labels_masked_YPARs.cmap		320 kbp	249 kbp	15.42	89.6%	186.663 X	7.7%	9.4%	Collected: 643.556 Gbp Target: 600.000 Gbp Percent: 107.259 %	643.556 Gbp
	3	Project: Chip QC Sample: 260701-0006153 Channel: Green 01 Enzyme: DLE-1 Reference: hg38_DLE1_0kb_labels_masked_YPARs.cmap		319 kbp	250 kbp	15.88	87.9%	178.007 X	7.9%	9.2%	Collected: 625.150 Gbp Target: 600.000 Gbp Percent: 104.192 %	625.150 Gbp

Figure 3. Run Metrics table.

- Plot:** The line color used in the plot for each flow cell. By deselecting the checkbox(s), a certain flowcell's data is hidden in the graph. Stratys chips have only one flowcell per chip.
- Flowcell:** The flowcell on the chip.
- Prep:** The names of the project, sample, channel, enzyme, and provided reference, as designated and entered by the user at Bionano Access.
- Avg N50>=150kbp:** The molecule N50 for all molecules that are ≥ 150 kbp in length. Given a set of molecules, the N50 is defined as the sequence length of the shortest molecule at 50% of the total molecule length.
- Avg N50>=20kbp:** The molecule N50 for all molecules that are ≥ 20 kbp in length.
- Avg Label Density per 100 kbp:** Average number of labels per 100 kbp for total DNA molecules that are ≥ 150 kbp in length.
- Avg Map Rate:** The average percentage of molecules (≥ 150 kbp in length) that map to the reference. If no reference is provided, the metric is 0%.
- Est Effective Coverage >=150 kbp:** The estimated effective coverage is calculated as follows: Avg Map Rate * Total DNA throughput/ length of the provided reference.
- Avg PLV (Average Positive Label Variance):** Percentage of molecule labels absent in the reference.
- Avg NLV (Average Negative Label Variance):** Percentage of reference labels absent in molecules.
- DNA Collected >=150 kbp:** The total amount of DNA ≥ 150 kbp that is collected for each flowcell during the run.
- Scan Count:** "Received" number reflects the number of scans that are submitted by the Instrument Control Software (ICS) to Bionano Access. "Mapped" number reflects the number of scans that have been converted to molecule data and had metrics generated.
- Cumulative DNA >=150 kbp:** "Collected" reflects the total amount of DNA that is collected in the flowcell across all runs of the chip. "Target" reflects the target throughput. "Percent" is the percentage of total collected throughput against the target throughput.

Cohort Data

To see the raw cohort data behind the graphs click the Cohort Data icon in the header. The system will display all the raw data points in a table that can be exported to Excel.




Figure 3a. Cohort Data icon

bionano

Access

Home > Chips > Dashboard > Cohort Data Points

Bionano Development | Administrator

Chip Run: 2b178fb7-1e51-4a46-aae4-8cd63cbbfe5e

Export to Excel

Flowcell	Scan	Bank	Cohort	Channel	Quantity	# Mals	NSD	NSD-20k	Mid SNR	Mid Inters	Label Date	Bpp	Stretch %	Scaling SD	Reference	# GoodMaps	Map Rate %	FP / 100kb	FP Rate %	FN Rate %	Coverage	SE	SF	SR
1	1	1	1	1	1.735e+8	857	278.982	210.068	291.3759	5041.259	19.572132	468	89.0	0.000e+0	3.083e+9	209	45.7	1.5	7.5	11.4	0.03	0.17	0.10	0.02
1	1	1	2	1	1.679e+8	805	301.002	218.803	302.5521	5197.7095	19.737543	468	88.7	1.207e-2	3.083e+9	189	45.4	1.3	6.7	9.9	0.02	0.18	0.09	0.02
1	1	2	1	1	2.322e+8	862	280.583	211.787	311.94815	5401.1304	20.416174	468	88.6	0.000e+0	3.083e+9	280	45.2	1.3	6.8	9.1	0.03	0.17	0.09	0.02
1	1	2	2	1	2.914e+8	1.145	284.337	197.731	322.18994	5574.8794	21.040806	467	88.9	2.567e-2	3.083e+9	281	36.3	1.5	8.0	9.1	0.03	0.17	0.08	0.02
1	1	3	1	1	1.251e+8	443	327.917	231.257	268.8152	4561.8613	18.117008	0	0.0	0.000e+0	0.000e+0	0	0	0.0	0.0	0.0	0.00	0.00	0.00	0.00
1	1	3	2	1	1.213e+8	427	304.418	241.843	284.95667	4851.7134	18.964114	0	0.0	0.000e+0	0.000e+0	0	0	0.0	0.0	0.0	0.00	0.00	0.00	0.00
1	1	4	1	1	1.799e+8	830	314.285	239.370	295.92383	5098.58	19.789196	468	88.6	0.000e+0	3.083e+9	234	52	1.4	7.4	8.2	0.03	0.18	0.09	0.02
1	1	4	2	1	2.244e+8	796	305.309	234.317	308.2146	5330.641	20.265171	469	88.6	0.000e+0	3.083e+9	271	48.5	1.3	6.8	9.0	0.04	0.18	0.10	0.02
1	2	1	1	1	6.152e+8	21.595	304.147	238.841	114.65316	3672.7932	16.54264	470	88.3	0.000e+0	3.083e+9	882	88.2	1.9	8.8	9.0	0.00	0.17	0.10	0.01
1	2	1	2	1	6.965e+8	23.752	313.076	263.837	125.12468	3842.421	16.28282	473	87.8	2.368e-2	3.083e+9	867	86.7	1.8	9.1	8.6	1.96	0.17	0.09	0.01
1	2	2	1	1	7.471e+8	25.948	307.966	254.995	113.89196	3272.1924	16.197618	473	87.7	2.401e-2	3.083e+9	861	86.1	1.7	8.7	8.4	2.08	0.18	0.09	0.01
1	2	2	2	1	6.370e+8	22.434	302.370	245.624	115.78664	3139.9132	16.056776	472	88.0	8.208e-3	3.083e+9	900	90	1.5	7.9	9.2	1.86	0.17	0.10	0.01
1	2	3	1	1	7.363e+8	24.956	317.398	266.848	88.489975	2982.8457	15.910230	472	87.8	2.814e-2	3.083e+9	885	88.5	1.5	7.9	8.4	2.11	0.18	0.09	0.01
1	2	3	2	1	7.295e+8	24.487	321.690	271.324	75.26127	2853.583	16.07249	474	87.6	1.419e-2	3.083e+9	908	90.8	1.5	7.8	8.6	2.14	0.17	0.10	0.01
1	2	4	1	1	6.637e+8	22.132	325.783	281.273	88.954285	2983.7895	15.995227	474	87.6	0.000e+0	3.083e+9	906	90.6	1.6	8.0	8.6	1.95	0.17	0.10	0.01
1	2	4	2	1	5.570e+8	16.755	320.969	277.297	100.45904	3041.5984	15.988075	473	87.7	2.113e-2	3.083e+9	887	88.7	1.5	7.7	8.8	1.60	0.17	0.10	0.01
1	3	1	1	1	6.195e+8	22.327	292.904	221.880	100.07415	3941.8338	16.988043	470	88.3	5.170e-3	3.083e+9	857	85.7	1.9	10.1	8.6	1.72	0.17	0.10	0.01
1	3	1	2	1	7.413e+8	25.457	311.760	256.421	113.32422	4016.396	16.433357	473	87.8	5.609e-3	3.083e+9	880	88	1.8	9.2	8.6	2.11	0.18	0.10	0.01
1	3	2	1	1	7.755e+8	26.854	309.896	255.676	107.35426	3756.1594	16.433357	473	87.7	0.000e+0	3.083e+9	837	83.7	1.7	8.5	8.6	2.10	0.18	0.10	0.01
1	3	2	2	1	6.701e+8	23.294	306.813	247.181	107.7314	3502.4277	16.106887	471	88.1	0.000e+0	3.083e+9	857	85.7	1.6	8.3	8.8	1.86	0.17	0.10	0.01
1	3	3	1	1	7.829e+8	26.907	311.625	255.174	88.2089	3211.9766	16.189346	472	87.8	1.609e-2	3.083e+9	871	87.1	1.6	8.2	8.4	2.21	0.17	0.10	0.01
1	3	3	2	1	7.989e+8	27.209	315.250	258.779	75.68981	3484.8713	16.182284	473	87.7	7.370e-3	3.083e+9	867	86.7	1.6	8.3	8.5	2.24	0.17	0.10	0.01
1	3	4	1	1	7.678e+8	25.975	318.730	261.621	93.89646	3954.5654	16.285066	473	87.7	1.346e-2	3.083e+9	859	85.9	1.6	8.5	8.9	2.14	0.18	0.10	0.01
1	3	4	2	1	6.750e+8	23.916	315.680	258.007	106.77955	3629.9463	16.200238	473	87.7	1.541e-2	3.083e+9	824	82.4	1.5	7.9	9.4	1.80	0.17	0.10	0.01
1	4	1	1	1	6.357e+8	22.440	298.642	233.642	105.29049	3840.812	16.056441	470	88.4	1.160e-2	3.083e+9	848	84.8	1.7	8.9	8.7	1.75	0.17	0.10	0.01
1	4	1	2	1	7.620e+8	25.752	316.714	262.814	105.58786	3686.8485	16.361032	472	87.9	0.000e+0	3.083e+9	858	85.8	1.7	8.7	8.8	2.12	0.17	0.10	0.01
1	4	2	1	1	7.681e+8	26.315	311.283	254.952	98.16874	3472.8716	16.25947	473	87.8	0.000e+0	3.083e+9	852	85.2	1.6	8.4	9.3	2.11	0.17	0.10	0.01
1	4	2	2	1	6.915e+8	22.844	308.650	246.567	97.26126	3382.5417	16.165916	471	88.1	5.220e-3	3.083e+9	855	85.5	1.6	8.4	9.6	1.83	0.17	0.10	0.01
1	4	3	1	1	7.726e+8	26.839	311.868	250.764	79.52334	3284.1028	16.09504	472	87.8	1.622e-2	3.083e+9	860	86	1.5	8.0	8.4	2.15	0.17	0.10	0.01
1	4	3	2	1	7.950e+8	27.896	316.882	254.205	73.13844	3548.2368	16.167381	473	87.7	0.000e+0	3.083e+9	871	87.1	1.6	8.3	8.8	2.24	0.18	0.10	0.01
1	4	4	1	1	7.803e+8	25.719	319.891	257.824	91.04515	3887.2418	16.322714	473	87.7	1.710e-2	3.083e+9	845	84.5	1.6	8.2	9.0	2.08	0.17	0.10	0.01
1	4	4	2	1	6.802e+8	22.836	315.626	252.109	99.801634	3616.8525	16.300024	473	87.7	1.357e-2	3.083e+9	816	81.6	1.6	8.1	9.1	1.77	0.17	0.10	0.01

Figure 4b. Cohort Data icon

Run Troubleshooting

Problem	Possible Causes	Recommended Actions
1. No data is shown in the Bionano Access dashboard.	Configuration issue in either ICS or Access.	Contact Bionano Technical Support (support@bionanogenomics.com).
2. No or low DNA throughput for molecules ≥ 150 kbp, starting from scan 1 to scan 5 (e.g., <10 Gbp per scan for human samples for Saphyr P/N 60239, <20 Gbp per scan for human samples for Saphyr P/N 60325 or Stratys P/N 60534).	Labeled DNA is <4 ng/ μ l (Direct Label and Stain - DLS) or <3 ng/ μ l (Nick, Label, Repair, and Stain - NLRS).	Repeat quantitation of labeled sample. Ensure labeled DNA concentration is within range: 4-16 ng/ μ l (DLS), 3-10 ng/ μ l (NLRS) (CV<0.30).
	Low N50.	See Row 10 of this table.
3. Decrease in throughput over time.	Chip clogging.	Check for labeled DNA homogeneity (CV<0.30). If DNA is not homogeneous, keep the DNA at room temperature overnight. Before loading on a new flowcell, use a regular p200 tip to pipette the labeled DNA up and down 2-3 times.
4. More than 20% decrease in map rate over time.	DNA sticks to the chip channels.	Check for labeled DNA homogeneity (CV<0.30). If DNA is not homogeneous, keep the DNA at room temperature overnight. Before loading on a new flowcell, use a regular p200 tip to pipette the labeled DNA up and down 2-3 times.
5. Map rate is zero or consistently low (<40%).	No reference or wrong reference provided.	Wait for the run to complete. After the Molecule object is imported to Access, edit the object so the correct reference is used.
	Low quality reference.	If possible, provide a better (more complete and/or more accurate) reference after run completion (see above).
	Lower than expected label density.	See Row 6 of this table.
	Wrong color selected.	Contact Bionano Technical Support (support@bionanogenomics.com).

Problem	Possible Causes	Recommended Actions
	Low throughput (<1 Gbp/scan).	See Row 2 of this table.
6. Average label density is lower than the expected label density (<i>in-silico</i> digestion). Greater than 3 labels/100 kbp difference between observed label density and expected label density would be considered out of range.	Inhibitory substances in the raw gDNA.	Increase Proteinase K digestion during DNA extraction and labeling, if applicable.
	Incorrect enzyme to DNA ratio in the labeling reaction (i.e., too much DNA and/or too little enzyme).	Ensure correct volumes of both DNA and enzyme are added to the labeling reaction.
	Others for DLE-1 labeled samples.	Please refer to Troubleshooting section of the <i>Bionano Prep Direct Label and Stain Protocol</i> (CG-30206).
7. Average label density is higher than the expected label density. Greater than 3 labels/100 kbp difference between observed label density and expected label density would be considered out of range.	Incorrect enzyme to DNA ratio in the labeling reaction (i.e., too little DNA and/or too much enzyme).	Ensure correct volumes of both DNA and enzyme are added to the labeling reaction.
8. High average NLV. NLV >20% (NLRS) or NLV >15% (DLS). *Only evaluate NLV when the map rate is consistent throughout the run and the average map rate is >40%.	Low enzyme activity, or inhibitory substances in the raw gDNA.	Check expiration of labeling enzyme and consider labeling a known control to evaluate enzyme labeling efficiency. Increase Proteinase K digestion during DNA extraction and labeling, if applicable.
9. High average PLV. PLV >15% (NLRS) or PLV >10% (DLS). *Only evaluate PLV when the map rate is consistent throughout the run and the average map rate is >40%.	Non-specific labels in the sample.	Repeat labeling and membrane adsorption. Wet the underside of DLS membrane with 1X DLE-1 buffer up to 10 minutes before sample application. Seal wells to prevent evaporation. Follow recommended incubation times.
10. Fragmented DNA (short molecules). Filtered (>150 kbp) N50 is < 230 kbp, or unfiltered (>20 kbp) N50 is < 130 kbp.	Poor quality of starting material.	Refer to Bionano released DNA extraction protocols for guidance.
	Improper handling of purified gDNA.	Avoid vortexing, rapid pipetting, or excessive pipetting with standard bore tips and use commercial wide-bore tips when appropriate.

Chip Metrics Graphs

Bionano Access enables users to view trends in chip run performance over time, from the **Chips** module under the **Metrics** tab. This view is useful for evaluating run consistency across many samples and flowcells over several weeks to months.

The chip metrics consist of six graphs (**Figure 4**): Total DNA (>=150 kbp), N50 (>=150 kbp), Average label density (>=150 kbp, Map rate (%), DNA per scan (Gbp), and Longest molecule (kbp). The x-axis for each graph is the date of data collection, and users may toggle to view a run performance period of 30, 60, or 90 days. By clicking on a graph, a user can zoom in, and then click again to zoom out.

Each datapoint represents one flowcell run on a Saphyr or Stratys instrument. Hovering over a given datapoint will display a tooltip of the corresponding y-axis value, run ID, and flowcell. Clicking on a data point will bring the user to the dashboard corresponding to that flowcell.




Figure 5. Chip metrics graphs

Technical Assistance

For technical assistance, contact Bionano Technical Support.

You can retrieve documentation on Bionano products, SDS's, certificates of analysis, frequently asked questions, and other related documents from the Support website or by request through e-mail and telephone.

TYPE	CONTACT
Email	support@bionano.com
Phone	Hours of Operation: Monday through Friday, 9:00 a.m. to 5:00 p.m., PST US: +1 (858) 888-7663 Monday through Friday, 9:00 a.m. to 5:00 p.m., CET UK: +44 115 654 8660 France: +33 5 37 10 00 77 Belgium: +32 10 39 71 00
Website	www.bionano.com/support
Address	Bionano , Inc. 9540 Towne Centre Drive, Suite 100 San Diego, CA 92121

Legal Notice

For Research Use Only. Not for use in diagnostic procedures.

This material is protected by United States Copyright Law and International Treaties. Unauthorized use of this material is prohibited. No part of the publication may be copied, reproduced, distributed, translated, reverse-engineered or transmitted in any form or by any media, or by any means, whether now known or unknown, without the express prior permission in writing from Bionano Genomics, Inc. Copying, under the law, includes translating into another language or format. The technical data contained herein is intended for ultimate destinations permitted by U.S. law. Diversion contrary to U. S. law prohibited. This publication represents the latest information available at the time of release. Due to continuous efforts to improve the product, technical changes may occur that are not reflected in this document. Bionano Genomics, Inc. reserves the right to make changes in specifications and other information contained in this publication at any time and without prior notice. Please contact Bionano Genomics, Inc. Customer Support for the latest information.

BIONANO GENOMICS, INC. DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. TO THE FULLEST EXTENT ALLOWED BY LAW, IN NO EVENT SHALL BIONANO GENOMICS, INC. BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT BIONANO GENOMICS, INC. IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

Patents

Products of Bionano Genomics® may be covered by one or more U.S. or foreign patents.

Trademarks

The Bionano logo and names of Bionano products or services are registered trademarks or trademarks owned by Bionano Genomics, Inc. ("Bionano") in the United States and certain other countries.

Bionano™, Bionano Genomics®, Saphyr®, Bionano Access™, Stratys™, and Stratys™ Chip are trademarks of Bionano Genomics, Inc. All other trademarks are the sole property of their respective owners.

No license to use any trademarks of Bionano is given or implied. Users are not permitted to use these trademarks without the prior written consent of Bionano. The use of these trademarks or any other materials, except as permitted herein, is expressly prohibited and may be in violation of federal or other applicable laws.

© Copyright 2026 Bionano Genomics, Inc. All rights reserved