



Bionano Solve™ v3.8.3.2 Release Notes

DOCUMENT NUMBER:

RNOTE-00007

DOCUMENT REVISION:

I

Effective Date:

Aug-27-2026

Table of Contents

Revision History	4
Introduction	5
Compatibility	5
References	5
Improvements for Solve 3.8	5
Support for New Reference Genomes	5
Singularity	5
Improved Support for Standard Genomics File Formats	6
Control Database	6
Gene Annotations	6
Y-PAR Masked Human References	6
Stable Region Analysis	7
Known Issues and Limitations	7
Improvements for Solve 3.8.1	7
Guided Assembly	7
RedHat 9 and Slurm Support	7
Additional Issues	7
Improvements for Solve 3.8.2	8
Quality Score Standardization	8
Copy Number Variant Quality Assessment	8
Other Improvements	8
Improvements for Solve 3.8.2.1	8

Compatibility for Gen 5 Saphyr compute servers	8
EnFocus Mode	9
Hybrid Scaffold Restored	9
BAM Conversion of Molecule to Reference Alignment Restored	9
Improvements for Solve 3.8.3	9
Restored SV Sensitivity	9
Restored SV Specificity	9
Increased Control Database	9
Updated SV and CNV ISCN notation with the 2024 version	10
Other Improvements	10
Improvements for Solve 3.8.3.1	10
Improvements for Solve 3.8.3.1Fixed missing dependency	10
Fixed missing tandem duplication at <i>KMT2A</i>	10
Improvements for Solve 3.8.3.2	11
Support for new compute hardware	11
Other improvements	11
Technical Assistance	11
Legal Notice	12
Patents	12
Trademarks	12

Revision History

REVISION	NOTES
A	Initial release of document.
B	Updated for Solve 3.8.1 release including information on Guided Assembly, RedHat9/Slurm support and issues addressed since Solve 3.8
C	Updated to add information under Known Issues and Limitations section.
D	Updated for Solve 3.8.2 release
E	Updated for Solve 3.8.2.1 update
F	Documenting additional fixes incorporated into Solve 3.8.2.1
G	Updated for Solve 3.8.3 release
H	Updated for Solve 3.8.3.1 hotfix
I	Updated for Solve 3.8.3.2 hotfix

Introduction

This document describes the release of Bionano Solve™ 3.8. This is an overview of the fixes and improvements of the Bionano Solve analysis tools and pipelines to provide a better understanding of the impact of moving to this version of the software. Should any questions arise, please contact support@bionano.com.

Bionano Tools and Bionano Solve are combined and branded as Bionano Solve. Bionano Solve is installed on Saphyr® Compute, Bionano Compute, and Bionano Access™ Servers before server shipment and installation.

Bionano Solve (folder “tools”) is located at the /home/bionano directory on the Compute server. The folder contains a collection of tools and scripts. Each individual tool is versioned independently. These tools together perform bioinformatics analyses on the Compute server.

Compatibility

Bionano Solve 3.8 is compatible with Bionano Access 1.8 only.

References

Visit <https://bionano.com/software-and-data-analysis-support-materials//> for file format specifications and Theory of Operation documents.

Improvements for Solve 3.8

Support for New Reference Genomes

TELOMERE-TO-TELOMERE CONSORTIUM REFERENCE

- Implemented system level support for the T2T-CHM13v2.0 reference genome
- Integrated T2T-CHM13v2.0 gene annotation
- Created Bionano Control Database versions for *de novo* Assembly and Rare Variant Analysis pipeline

MM39 REFERENCE

- Implemented system level support for the mouse mm39/GRCm39 reference genome
- Integrated mm39 gene annotation
- Created Bionano Control Database versions for *de novo* Assembly and Rare Variant Analysis pipeline

Singularity

- Adopted Singularity for dependency management of the Solve pipeline (replacing Docker)
- Provides method for portable, consolidated dependency management across supported Compute hardware

- Designed for ease-of-use on high-performance computing (HPC) systems
- Provides improved security for customers by removing requirements for running with elevated or root user privileges

Improved Support for Standard Genomics File Formats

- Introduced OGM BAM output (binary version of the Sequence Alignment/Map format) for molecule-to-reference alignments.
- Updated OGM Variant Call File (VCF) with improved representations of variants and quality filters.
- Updated OGM VCF to include Solve AOH calls for *de novo* and Guided Assembly – constitutional analyses.
- Updated OGM VCF to contain clustered variant calls only to remove redundant calls.
 - Translocations not included in clustering for initial release.
- Enables import of OGM data into VIA as well as any standards-compliant genomics software that supports BAM and VCF.

Control Database

- Added DLE-1 samples to human control database for a total of 285 individuals.
- Complete reanalysis of control databases for *de novo* Assembly and Rare Variant Analysis for hg19 and hg38.
- Introduction of control database for T2T-CHM13v2 reference genome.

Gene Annotations

- Updated gene annotation used by Variant Annotation Pipeline to harmonize gene annotations used by Access and VIA.
- Updated gene annotation to latest RefSeq builds for hg38 and hg19.
- Introduced gene annotation for T2T-CHM13v2.0 and mm39.
- Streamlined gene annotation to provide a single, high quality, current annotation for human analyses while allowing use of custom annotations if desired.

Y-PAR Masked Human References

- Added versions of human reference genomes that mask out the pseudoautosomal regions (PAR) on chromosome Y as do approaches used in many Next-Generation Sequencing (NGS) analysis workflows.
- This is done to address the sequence homology in these regions with the corresponding regions on chromosome X which can interfere with map and molecule alignments to the reference. Masking of these regions has been shown to improve structural and copy number variant calling for genes such as *CRLF2* that are in or near the region.
- Masked references are provided as options for hg19 and hg38 and by default for T2T-CHM13v2.0. For best results with SV detection performance and annotation with the Bionano control databases, the masked references are recommended for hg19 and hg38.

Stable Region Analysis

- Incorporated analysis of stable regions as a quality control step.
- Initially implemented in the EnFocus™ FSHD and EnFocus™ Fragile X applications, assessment of stable regions is now included in the *de novo* and Guided Assembly whole genome analysis pipelines.
- Output included in informatics report

Known Issues and Limitations

- Control databases were generated using the Y-PAR masked versions of hg19 and hg38. These databases can be used to annotate data analyzed using the unmasked versions of the reference genome, but SV calls in Y PAR1 and PAR2 will show as having 0% match in the control database.

Improvements for Solve 3.8.1

Guided Assembly

- Introduction of Guided Assembly pipeline as an alternative to *de novo* assembly and Rare Variant Analysis
- Mode of *de novo* Assembly that uses the reference genome as seed for refinement and extension
- Constitutional analysis targeted for 400 Gb and 800 Gb coverage tiers
- Low Allele Fraction (LAF) – targeted for 1.5 Tb coverage tier for cancer analysis
- Improves variant detection in low-allele fraction samples with improvements in insertions/deletions and duplication detection
- Dedicated control database for both modes for hg19, hg38 and T2T-CHM13v2.0 reference genomes

RedHat 9 and Slurm Support

- Established support for RedHat 9 Linux operating system. RedHat Linux is a commercially supported distribution and is offered as a replacement to CentOS Linux 7 which will reach end of life (EOL) on June 30, 2024.
- Implemented use of the Slurm Workload Manager for all analysis pipelines. Slurm replaces Sun Grid Engine in the RedHat 9 version of the Solve pipeline.
- CentOS and SGE installers are provided for Solve 3.8.1.

Additional Issues

- Corrected defect in dual and trio modes of Variant Annotation Pipeline where molecule support was not calculated for case or parental samples.
- Corrected issue where CNV calls for mouse reference mm10 reported with confidence scores of zero.
- Corrected issue so that translocation calls are clustered in the VCF as are other variant types.
- Updated representation of translocations in VCF to reflect directionality of breakpoints more accurately.
- Corrected confidence intervals for translocations in VCF so that each breakpoint is assessed separately.

- Remove redundant samples from human DLE-1 control databases. Solve 3.8 reported 384 individuals while number of unique individuals is 285.
- Corrected Variant Annotation Pipeline human gene annotation files to contain all transcripts rather than only primary transcript.

Improvements for Solve 3.8.2

Quality Score Standardization

- Implemented uniform quality scoring in VCF for structural variants. High confidence variants of any type are defined as variants with a QUAL score ≥ 20 .
- Implemented confidence score for insertion and deletions 300bp – 1kbp.
- Implemented confidence score for duplication structural variants.

Copy Number Variant Quality Assessment

- Added coverage quality assessment to detect presence of telomeric gain artifacts. Assessment reported in the informatics report of all human and mouse whole genome pipelines as well as EnFocus™ FSHD and EnFocus™ Fragile X.

Other Improvements

- Fixed issue with HybridScaffold perl locale warnings.
- Corrected VCF representation of deletion calls produced by the RVA pipeline to show genotype as undetermined rather than defaulting to heterozygous.
- Fixed bug in de novo assembly pairwise alignment stage that in some cases increased fragmentation of assembly. While structural variant calls do not appear affected, total number of genome maps was elevated in some cases.
- Refreshed de novo assembly control databases for hg19, hg38 and T2T-CHM13v2.0 reference genomes.
- Added improvements to pipeline reproducibility to reduce variation between analyses.
- Added improvements to pipeline robustness to remove silent failures and improve reporting of root cause errors.
- Updated handling of molecule sizing errors and molecules with outlying label density to improve compute performance.

Improvements for Solve 3.8.2.1

Compatibility for Gen 5 Saphyr compute servers

- Periodically we update the generation of our compute servers for the Saphyr system to leverage the latest technologies. This update adds cluster arguments and other small changes needed to support the Saphyr

Compute Server, Gen5 and Bionano Compute Server, Gen5. Use of the Gen5 compute servers requires running Solve 3.8.2.1 and the Red Hat 9 operating system.

EnFocus Mode

- This version has been updated to support running EnFocus FSHD and Fragile-X jobs running directly on the Bionano Access Servers that have SLURM installed.

Hybrid Scaffold Restored

- The Hybrid Scaffold operation in Solve 3.8.2 was not functional. This version restores operation of the Hybrid Scaffold analysis.

BAM Conversion of Molecule to Reference Alignment Restored

In Solve 3.8.2, the conversion step would fail in the rare case when there is no molecule aligned to a reference map. This version restores that operation.

Improvements for Solve 3.8.3

Restored SV Sensitivity

- Enhanced sensitivity of small insertions and deletions during SV detection. An error was caught in the computation of confidence score with more stringent sizing error parameters, and consequently, sensitivity is restored. In addition, by removing an antiquated duplicate filter, translocations with multiple match groups on the same map can now be detected.

Restored SV Specificity

- This version has been restored to remove low confidence insertions and deletions, < 500 bp, in the de novo and guided assemblies. Also, other false insertions and deletions < 1 kbp were removed by checking neighboring label intervals and fine-tuning scoring parameters. The specificity of translocation calls was improved by increasing the alignment P-Value from 1e-8 to 1e-10 in the rare variant analysis and guided assembly - LAF, and by suppressing calls with long simple repeats in the intervening region between breakpoints.

Increased Control Database

- The control database has been increased to 341 individuals, while SV mask regions are updated with the added samples.

Updated SV and CNV ISCN notation with the 2024 version

- All SV and CNV are annotated with the 2024 version of the International System for Human Cytogenomic Nomenclature (ISCN) standard. A bug in the order of chromosome number and band in translocation calls has been fixed.

Other Improvements

- Fixed rare variant analysis missing VAF regression
- Fixed the missing files in pipeline result zip file
- Fixed BAM conversion error when a reference chromosome has no molecule alignments
- Corrected Bonferroni adjustment of BAM MAPQ scores
- Fixed molecule misalignments in refine final and restored the correct zygoty label
- Added priority feature to SLURM configuration, pipeline and launch scripts for Stratys Compute
- Removed extraneous files from assembly result zip files to reduce assembly package size
- Added elapsed and CPU time to informatics report
- Added split alignment code to alignmolref to enable faster pipeline runtime (disabled by default)
- Included sex information in informatics report
- Updated gene and other feature files in Bionano Solve and VIA
- Bypass BAM conversion when a reference contains chromosomes larger than $2^{29}-1$ bp, the limit of BAI index
- Changed Solve version string to 3.8.3 in multiple pipelines

Improvements for Solve 3.8.3.1

Improvements for Solve 3.8.3.1Fixed missing dependency

- Our singularity container was missing the rlog package. The pipeline automatically will install the missing dependency during runtime, but it doesn't work at dark sites where the compute node does not have internet access. The singularity image was updated to include rlog.

Fixed missing tandem duplication at *KMT2A*

- Adjustments were made to detect missing tandem duplications at *KMT2A*. Tests were performed to confirm existing variant detection performance was not affected. Control databases were also updated to account for the adjustments made.

Improvements for Solve 3.8.3.2

Support for new compute hardware

- Added configurations to accelerate pipeline operations for the new computes

Other improvements

- Added routine check to remove GPU deadlocks. This is a known but very rare issue with Bionano Solve running on computes with GPUs, where the GPU lock used to control GPU access could not be released. That leads to stalled processes.
- Removed the calling of false terminal deletions in guided assembly - constitutional

Technical Assistance

For technical assistance, contact Bionano Technical Support.

You can retrieve documentation on Bionano products, SDS's, certificates of analysis, frequently asked questions, and other related documents from the Support website or by request through e-mail and telephone.

TYPE	CONTACT
Email	support@bionano.com
Phone	Hours of Operation: Monday through Friday, 9:00 a.m. to 5:00 p.m., PST US: +1 (858) 888-7663 Monday through Friday, 9:00 a.m. to 5:00 p.m., CET UK: +44 115 654 8660 (UK) France: +33 5 37 10 00 77 Belgium: +32 10 39 71 00
Website	www.bionano.com/support
Address	Bionano, Inc. 9540 Towne Centre Drive, Suite 100 San Diego, CA 92121

Legal Notice

For Research Use Only. Not for use in diagnostic procedures.

This material is protected by United States Copyright Law and International Treaties. Unauthorized use of this material is prohibited. No part of the publication may be copied, reproduced, distributed, translated, reverse-engineered or transmitted in any form or by any media, or by any means, whether now known or unknown, without the express prior permission in writing from Bionano Genomics, Inc. Copying, under the law, includes translating into another language or format. The technical data contained herein is intended for ultimate destinations permitted by U.S. law. Diversion contrary to U. S. law prohibited. This publication represents the latest information available at the time of release. Due to continuous efforts to improve the product, technical changes may occur that are not reflected in this document. Bionano Genomics, Inc. reserves the right to make changes to specifications and other information contained in this publication at any time and without prior notice. Please contact Bionano Genomics, Inc. Customer Support for the latest information.

BIONANO GENOMICS, INC. DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. TO THE FULLEST EXTENT ALLOWED BY LAW, IN NO EVENT SHALL BIONANO GENOMICS, INC. BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT BIONANO GENOMICS, INC. IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

Patents

Products of Bionano Genomics® may be covered by one or more U.S. or foreign patents.

Trademarks

The Bionano logo and names of Bionano products or services are registered trademarks or trademarks owned by Bionano Genomics, Inc. ("Bionano") in the United States and certain other countries.

Bionano™, Bionano Genomics®, Bionano Solve™, Saphyr®, Saphyr Chip®, Bionano Access™, VIA™ software, and Bionano EnFocus™ are trademarks of Bionano Genomics, Inc. All other trademarks are the sole property of their respective owners.

No license to use any trademarks of Bionano is given or implied. Users are not permitted to use these trademarks without the prior written consent of Bionano. The use of these trademarks or any other materials, except as permitted herein, is expressly prohibited and may be in violation of federal or other applicable laws.

© Copyright 2026 Bionano Genomics, Inc. All rights reserved.