



# **Bionano VIA™ Software Version 7.2.2 Release Notes**

**DOCUMENT NUMBER:**

RNOTE-00016

**DOCUMENT REVISION:**

D

**Effective Date:**

Aug-27-2026

---

## Table of Contents

|                                    |           |
|------------------------------------|-----------|
| <b>Revision History</b>            | <b>3</b>  |
| <b>Bionano VIATM Software</b>      | <b>4</b>  |
| Introduction                       | 4         |
| Supported Implementation Options   | 4         |
| Compatibility                      | 5         |
| <b>New Features in Version 7.2</b> | <b>5</b>  |
| Specific Region Sample Search      | 7         |
| Tickets                            | 10        |
| <b>Technical Assistance</b>        | <b>16</b> |
| <b>Legal Notice</b>                | <b>17</b> |
| Patents                            | 17        |
| Trademarks                         | 17        |

## Revision History

| REVISION | NOTES                        |
|----------|------------------------------|
| A        | Initial release.             |
| B        | Updates for VIA 7.2 release  |
| C        | Updates for VIA 7.2.1 update |
| D        | Updates for VIA 7.2.2 update |

# Bionano VIA™ Software

This document describes the v7.2.1 release of VIA software. This document provides an overview of what is changing with this release so that users may better understand the impact of moving to this version. Should users have any questions please contact [software-support@bionano.com](mailto:software-support@bionano.com).

## Introduction


Variant Intelligence Applications™ (VIA) software is a complete and integrated solution for the visualization, interpretation and reporting of genomic variants from multiple technology types. By supporting multiple genome-wide data modalities, VIA provides the most comprehensive view of genomic variants of any interpretation, annotation, and reporting software tool available. As a platform-agnostic tertiary analysis solution, VIA stores and manages distinct types of genomic data from various platforms enabling the extraction of meaningful insights from standalone or combined analysis.

The software includes algorithms to detect copy number variants (CNV) from major microarray vendors, optical genome mapping (OGM), and next generation sequencing (NGS) methodologies as well as Absence of Heterozygosity (AOH), from data types that assess B-allele frequency. VIA also provides interpretation assistance to analyze CNVs, Loss of Heterozygosity (LOH) and Structural Variants (SV) from OGM data. As a centralized analysis solution spanning technologies and application areas, VIA provides an efficient environment to keep pace with advancements in technology while retaining access to historical platform data. By being adaptive to whichever technology is used to generate CNV, LOH, or SV genomic variants, VIA provides rich annotations for the co-analysis of sequence variants from NGS to provide a complete picture of genomic variation and reveal more answers for disease association.

## Supported Implementation Options

VIA is a scalable enterprise solution comprised of three components: the VIA Server, VIA Processing, and VIA Client applications. The server is the primary component of the system that manages the data in the repository. Installed adjacent to the Server is the Processing application that executes processing jobs for segmentation and annotation. The Client software is the User Interface (UI) utilized by all users to access sample data hosted in the VIA Server. There are three implementation options to deploy the VIA system; installation of the VIA Server and VIA Processing on the Bionano Access Server (BAS), local network server, or hosted in a cloud environment as illustrated below in **Figure 1**. In each implementation option, the VIA Client software is installed locally on each analyst's workstation with network access to the VIA Server. Further information on hardware requirements is available in the *Bionano VIA System Requirements* guides (CG-30577 and CG-30580).

- Option 1: VIA Server and VIA Processing are hosted on the BAS and is a common installation strategy for OGM users primarily interested in analyzing OGM data in VIA software.
- Option 2: VIA Server and VIA Processing are hosted on a local network and is a common installation strategy for NxClinical users primarily interested in array or NGS technology.
- Option 3: VIA Server and VIA Processing are hosted on a site-administered cloud environment and is a common installation strategy for users seeking production-scale analysis with VIA.



**Figure 1.** VIA implementation options illustration

## Compatibility

VIA 7.2.1 is compatible with Bionano Access™ 1.8.3 and the Bionano Solve® v3.8.3 pipeline running on either Saphyr® Compute, Bionano Compute servers or Bionano Compute On Demand v1.4.3. Please refer to the Bionano VIA System Requirements guide (CG-30577) for hardware requirements to install the software on local network and client workstations.


## New Features in Version 7.2

### Similar Previous Cases for Structural Variants

VIA now supports similar previous cases for structural variants. In the track view if similar previous cases are available a carrot (triangle) will be displayed just the left variant type. If you click on the carrot, the system will display similar cases as shown below. Click the carrot again to collapse the similar previous cases. You can click on a similar previous case name in the track view to open that sample.



Similar previous cases for structural variants also appear in the Variant Details (image below), Event Table, and SV Filtering. Similar previous cases for decision tree classification of structural variants will be supported in future versions.



### Specific Region Sample Search

On the home page users can search for samples containing structural variants at specific regions or that overlap specific genes. Here are some examples:

*Insertion:EGFR*

*Insertion:chr7:100000-150000*

The structural variant types that are supported include the following:

- unpaired\_inversion
- paired\_inversion
- interchr\_translocation
- intrachr\_fusion
- inverted\_duplication
- sv-deletion
- sv-insertion
- sv-tandem\_duplication
- sv-inversion
- inversion\_breakpoint

## SAP Score for SV Events

The Significance Associated with Phenotype (SAP) score is now calculated for structural variant event types. The SAP Score column is available in the Event Table and can be used in Decision Trees to classify variants. For more information about SAP scores you can refer to this article:

<https://www.sciencedirect.com/science/article/abs/pii/S1525157824000291>

| Table Deleted Events Whole Genome Aneusomy Report |                  |                                     |           |                        |         |
|---------------------------------------------------|------------------|-------------------------------------|-----------|------------------------|---------|
| Select                                            | Event            | Chromosome Region                   | SAP Score | Similar Previous Cases | Cyt     |
| <input type="checkbox"/>                          | Loss             | chr2:61,055,001-62,076,000          | 1         | 0.97% (7/718)          | 2p11.2  |
| <input type="checkbox"/>                          | Gain             | chr2:89,450,001-90,369,000          | 1         | 2.37% (17/718)         | 2p11.2  |
| <input type="checkbox"/>                          | Gain             | chr10:47,434,212-48,090,211         | 1         | 5.01% (36/718)         | 10q24.3 |
| <input type="checkbox"/>                          | Deletion         | chr1:223,553,312-223,615,819        | 1         | 51.81% (372/718)       | 1q44    |
| <input type="checkbox"/>                          | Deletion         | chr4:99,083,117-99,129,639          | 1         | 0.97% (7/718)          | 4q21    |
| <input type="checkbox"/>                          | Insertion        | chr4:190,016,392-190,030,545        | 1         | 0.97% (7/718)          | 4q31.2  |
| <input type="checkbox"/>                          | Deletion         | chr10:47,776,804-47,870,426         | 1         | 46.24% (332/718)       | 10q24.3 |
| <input type="checkbox"/>                          | Deletion         | chr11:70,950,863-71,086,672         | 1         | 50.56% (363/718)       | 11q23   |
| <input type="checkbox"/>                          | Deletion         | chr11:87,972,433-88,007,271         | 1         | 55.15% (396/718)       | 11q23   |
| <input type="checkbox"/>                          | Deletion         | chr11:117,312,178-117,342,581       | 1         | 0.97% (7/718)          | 11q23   |
| <input type="checkbox"/>                          | Insertion        | chr12:20,319,274-20,322,819         | 1         | 0.97% (7/718)          | 12p11.2 |
| <input type="checkbox"/>                          | Deletion         | chr12:25,493,950-25,503,013         | 1         | 0.97% (7/718)          | 12p11.2 |
| <input type="checkbox"/>                          | Deletion         | chr13:113,652,332-113,753,146       | 1         | 51.81% (372/718)       | 13q31.2 |
| <input type="checkbox"/>                          | Deletion         | chr17:443,876-520,107               | 1         | 51.67% (371/718)       | 17p11.2 |
| <input type="checkbox"/>                          | Interchr Tran... | chr3:197,762,504-197,774,812;chr... | 1         | 0.97% (7/718)          | 3q21.31 |
| <input type="checkbox"/>                          | Tandem Dupl...   | chr17:2,250,783-2,389,481           | 1         | 0.97% (7/718)          | 17p11.2 |
| <input type="checkbox"/>                          | Interchr Tran... | chr3:197,604,666-197,630,064;chr... | 1         | 0.97% (7/718)          | 3q21.31 |
| <input type="checkbox"/>                          | Deletion         | chr19:14,583,236-14,596,858         | 1         | 0.97% (7/718)          | 19p11.2 |
| <input type="checkbox"/>                          | Tandem Dupl...   | chr21:10,295,808-10,468,228         | 1         | 0.97% (7/718)          | 21p11.2 |

| Sample Attributes |                                                                | Workflow | Event Classification | Decision Trees | Sample Review Preferences   |
|-------------------|----------------------------------------------------------------|----------|----------------------|----------------|-----------------------------|
| +                 |                                                                | -        | ✎                    | sap            | Apply (1 processed samples) |
| 1                 | CASE {ANY_CN_EVENT_KIND} {                                     |          |                      |                |                             |
| 2                 | IF { SCORE(SAP) < 0.12 } THEN {CLASSIFY("Likely Pathogenic")}} |          |                      |                |                             |



### Add ISCN for Normal Samples

In this release it is possible to add ISCN values to normal samples that have no events to report. The ISCN values generated would match the table below.

|        | <b>array</b>                | <b>ogm</b>                  |
|--------|-----------------------------|-----------------------------|
| female | <u>arr</u> (X,1-22)x2       | <u>ogm</u> (X,1-22)x2       |
| male   | <u>arr</u> (X,Y)x1,(1-22)x2 | <u>ogm</u> (X,Y)x1,(1-22)x2 |

Affected Status: ▼

ISCN:

☒ Add Normal Ploidy ISCN

### New SV Quality Filtering

In this release it is possible to filter structural variants based on their quality value as show below.

Filter Parameters ×

**SV QUALITY**


☐ Remove events with SV quality less than

### Support for Illumina Connected Annotation Files

Version 7.2 supports the import of data in Illumina Connected Annotations JSON format. For help creating a sample type to support this data format please contact our support team.

Mutli-Sample View Improvement

The mutli-sample view now sports an extra tab to show the ‘CN Probe Moving Average’ (image below). This new tab displays the average CN value for each sample along the genome. The left column includes a legend so you can determine which color matches to which sample.



Tickets

The features released with version 7.2 are listed in **Table 1**.

**Table 1.** Released features in v7.2

| Summary                                                                                    | Ticket  |
|--------------------------------------------------------------------------------------------|---------|
| Show the matched sample count of previous similar cases on the SV track                    | NC-5804 |
| Show similar previous cases in the table                                                   | NC-5946 |
| API to retrieve sample types from Access                                                   | NC-6216 |
| Enable tags for QC metrics so that these can be added and populated in the report template | NC-6311 |
| Support reporting for genes from databases used to annotate events                         | NC-6312 |
| Add the OGM chip serial number as standard attribute for OGM Sample Type                   | NC-6334 |

|                                                                                                |         |
|------------------------------------------------------------------------------------------------|---------|
| Display Similar previous cases dial for SV events                                              | NC-6343 |
| Modify SV index to store/use different overlap intervals based on event type                   | NC-6355 |
| Connect the REST API to the SV index via QueryManager to retrieve number of matched SV samples | NC-6358 |
| Add new SV previous cases parameters to the sample query menu                                  | NC-6365 |
| Allow export of HGVS nomenclature for SeqVars in the «T:Whole Genome Results» table            | NC-6370 |
| Add similar to previous cases filter for SV events                                             | NC-6396 |
| Connect the REST API to the SV index to retrieve samples with similar case data                | NC-6397 |
| Add new sort-order option "Sample Name" for SV previous similar cases                          | NC-6400 |
| Add query breakend direction to inputs for gene disruption queries parameterization            | NC-6401 |
| Display a CN probe moving average line for a set of samples in one chart                       | NC-6405 |
| Support export and copy to clipboard for the multisample view image                            | NC-6417 |
| Allow VIA to fall back to a default CRAM reference name when the UR field is missing.          | NC-6426 |
| Implement new DT Function SV_MOLECULE_COUNT                                                    | NC-6472 |
| Implement new DT Function SV_VAF                                                               | NC-6473 |
| Chromosome Region Serialize/Deserialize refactor for 7.2                                       | NC-6481 |
| Calculate a SAP score for SV event types                                                       | NC-6487 |
| export SV data in JSON                                                                         | NC-6493 |
| Create methods to access the SV index                                                          | NC-6517 |
| Vertical zoom for the multi sample CN moving average                                           | NC-6520 |
| Ignore corrupted samples from the homepage                                                     | NC-6534 |
| Allow users to search for samples containing SVs at specific regions or gene from the homepage | NC-6538 |
| Render expanded similar cases events in SV track                                               | NC-6542 |
| Allow users to include custom region values for an event in the word report                    | NC-6545 |
| Display the sample name in the CN prove moving average tab only if there are probes available  | NC-6548 |

|                                                                                                                                 |         |
|---------------------------------------------------------------------------------------------------------------------------------|---------|
| Allow sv track query sample reference to display various attributes                                                             | NC-6560 |
| Create new criteria using new header and decouple Gson parser version from Illumina Schema version.                             | NC-6564 |
| Disambiguate DbSnp parsing in ICA                                                                                               | NC-6565 |
| Parse transcript format flat map with source field in ICA                                                                       | NC-6566 |
| Only show SV prev sim cases query tab for samples with SVs                                                                      | NC-6569 |
| Perform aneusomy calculation for Affy CyChp data type                                                                           | NC-6582 |
| SV Previous Similar Cases Tool Tip                                                                                              | NC-6589 |
| Display classification colored bar on SV track and Similar Previous Cases                                                       | NC-6590 |
| Double right click on the probe moving average line plot should reset the plot                                                  | NC-6594 |
| Add a Y-axis to the probe moving average line plot of the multi-sample view                                                     | NC-6595 |
| Sv Previous Similar Cases options in admin preferences menu                                                                     | NC-6599 |
| SV previous similar case query should use confidence interval                                                                   | NC-6605 |
| Add ability to specify normal ISCN for samples regardless of test type                                                          | NC-6608 |
| Allow user to click on the label of a similar previous case to open the sample in a new tab                                     | NC-6617 |
| Change Clinvar.Significance into a string array instead of single string                                                        | NC-6624 |
| Remove ICA data version check (but keep the step name and the filter options)                                                   | NC-6625 |
| Clinvar VCV links should link out to variant page                                                                               | NC-6628 |
| Create a set of similar previous cases setting for Deletion, Duplication and another set for Insertion and Inverted duplication | NC-6629 |
| Add percent overlap parameter to similar previous case setting for Deletion and Duplication                                     | NC-6630 |
| Avoid unnecessary listener notifications for the SV similar cases preferences menu                                              | NC-6636 |
| SV track mouse for closed track plus text label font size                                                                       | NC-6637 |
| Graceful start on Ignite connection failure                                                                                     | NC-6696 |

The defects that have been addressed in Bionano Access version 1.8 are listed in **Table 2**.

**Table 2.** Defects Addressed in v7.2

| Summary                                                                                                                                                 | Ticket  |
|---------------------------------------------------------------------------------------------------------------------------------------------------------|---------|
| Compile DT Scripts in UTF-8 encoding                                                                                                                    | NC-6149 |
| Server reloads license without checking request key                                                                                                     | NC-6392 |
| When only using a VCF or only concerned with SeqVar, the list of the genes in panel will not appear on export                                           | NC-6467 |
| Prevent admin from adding attribute that contains only whitespace                                                                                       | NC-6512 |
| SCORE(EVIDENCE) for SV events should be prevented with useful message                                                                                   | NC-6514 |
| OGM BAF sample clustering exception when no probes at all                                                                                               | NC-6524 |
| Gene annotations should be cleared when a sample is duplicated or reset                                                                                 | NC-6531 |
| SeqVar processing types should not be selectable if the seqvar modality is not selected for the sample type                                             | NC-6541 |
| Attempting to navigate to the sample type associated with a sample from the homepage generates an error                                                 | NC-6550 |
| Parent of origin probe coloring should show up in whole genome view tab initially at sample opening                                                     | NC-6552 |
| Completion Exception - undefined manufacturer generated when adding a sample type with sample classes Array only, Methylation, NGS & Array and GxA-Cyto | NC-6555 |
| Address secondary matching edge cases                                                                                                                   | NC-6557 |
| Table columns such as parent of origin, inheritance comment fail to update if an event is manually edited or created                                    | NC-6558 |
| Memory leak in dpiScalingPreferences listeners                                                                                                          | NC-6559 |
| Unable to login due to error in sample type preference migration                                                                                        | NC-6561 |
| Sample order in the multi-sample view window should follow the sample order on home page                                                                | NC-6572 |
| Synchronize column width in main table and preferences table layout                                                                                     | NC-6584 |
| Failure to query samples from the homepage based on sequence variants at a particular location                                                          | NC-6596 |
| Error using «T:Whole Genome Results» tag with single-position Seqvar events                                                                             | NC-6604 |
| OGM label density systematic correction apply to all probes                                                                                             | NC-6611 |
| KB Relevant Genes Dialog exception with blank notes                                                                                                     | NC-6612 |
| Admin decision tree columns re-calculation SeqVar index bounds check                                                                                    | NC-6615 |

|                                                                                                    |         |
|----------------------------------------------------------------------------------------------------|---------|
| Exception manually adding CNV events to SeqVar samples                                             | NC-6616 |
| Seqvar samples not processing with VCF proc type                                                   | NC-6631 |
| Decision tree SeqVar event fix for ANY_SEQVAR_EVENT_KIND predicate                                 | NC-6638 |
| BAF Segment Value isn't being calculated for allelic events                                        | NC-6639 |
| cnScaledProbes and cnScaledAvg temp files should be deleted along with normal probe and snp files. | NC-6645 |
| Nirvana JSON loading accept "Illumina Annotation Engine" annotator                                 | NC-6653 |
| Column preference UI ArrayIndexOutOfBoundsException                                                | NC-6654 |
| Add Logging for Ignite                                                                             | NC-6661 |
| exception when loading older duo samples                                                           | NC-6667 |
| mouse wheel scrolls in the wrong direction                                                         | NC-6669 |
| SeqVar JSON omitted genotype field interpreted as dot (.)                                          | NC-6680 |
| ISCN correction from alpha testers                                                                 | NC-6681 |
| Similar cases score for multiple matches from same sample                                          | NC-6686 |
| Decision tree throwing error 22                                                                    | NC-6691 |
| Change Ignite Configuration                                                                        | NC-6695 |

**Table 2.** Defects Addressed in v7.2.1

| Summary                                    | Ticket  |
|--------------------------------------------|---------|
| Incorrect ISCN generated for Normal Ploidy | NC-6768 |
| Thread performance issue                   | NC-6756 |
| S3 batch loading .crai suffix fix          | NC-6792 |

**Table 3.** Defects Addressed in v7.2.2

| Summary | Ticket |
|---------|--------|
|---------|--------|

Latest ICS JSON file has formatting changes that impact import into VIA

NC-6844

## Known Issues

We have observed in some cases the VIA client is not allocated sufficient memory (NC-6871). We plan to update the memory settings in the VIA Client installer in the next release. For help adjusting the memory configuration for your existing VIA Client please contact support.

There is a known issue with IDAT Methylation probes on ChrX and ChrY (NC-6839). These issues will be resolved in the next release.

## Technical Assistance

For technical assistance, contact Bionano Technical Support.

You can retrieve documentation on Bionano products, SDS's, certificates of analysis, frequently asked questions, and other related documents from the Support website or by request through e-mail and telephone.

| TYPE    | CONTACT                                                                                                                                                                                                                                                        |
|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Email   | <a href="mailto:support@bionano.com">support@bionano.com</a> , <a href="mailto:software-support@bionano.com">software-support@bionano.com</a>                                                                                                                  |
| Phone   | <p>Hours of Operation:<br/> Monday through Friday, 9:00 a.m. to 5:00 p.m., PST<br/> US: +1 (858) 888-7663</p> <p>Monday through Friday, 9:00 a.m. to 5:00 p.m., CET<br/> UK: +44 115 654 8660<br/> France: +33 5 37 10 00 77<br/> Belgium: +32 10 39 71 00</p> |
| Website | <a href="http://www.bionano.com/support">www.bionano.com/support</a>                                                                                                                                                                                           |
| Address | <p>Bionano, Inc.<br/> 9540 Towne Centre Drive, Suite 100<br/> San Diego, CA 92121</p>                                                                                                                                                                          |



## Legal Notice

### **For Research Use Only. Not for use in diagnostic procedures.**

This material is protected by United States Copyright Law and International Treaties. Unauthorized use of this material is prohibited. No part of the publication may be copied, reproduced, distributed, translated, reverse-engineered or transmitted in any form or by any media, or by any means, whether now known or unknown, without the express prior permission in writing from Bionano Genomics, Inc. Copying, under the law, includes translating into another language or format. The technical data contained herein is intended for ultimate destinations permitted by U.S. law. Diversion contrary to U. S. law prohibited. This publication represents the latest information available at the time of release. Due to continuous efforts to improve the product, technical changes may occur that are not reflected in this document. Bionano Genomics, Inc. reserves the right to make changes in specifications and other information contained in this publication at any time and without prior notice. Please contact Bionano Genomics, Inc. Customer Support for the latest information.

BIONANO GENOMICS, INC. DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. TO THE FULLEST EXTENT ALLOWED BY LAW, IN NO EVENT SHALL BIONANO GENOMICS, INC. BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT BIONANO GENOMICS, INC. IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

### **Patents**

Products of Bionano Genomics® may be covered by one or more U.S. or foreign patents.

### **Trademarks**

The Bionano logo and names of Bionano products or services are registered trademarks or trademarks owned by Bionano Genomics, Inc. ("Bionano") in the United States and certain other countries.

Bionano™, Bionano Genomics®, Saphyr®, and VIA™ are trademarks of Bionano Genomics, Inc. All other trademarks are the sole property of their respective owners.

No license to use any trademarks of Bionano is given or implied. Users are not permitted to use these trademarks without the prior written consent of Bionano. The use of these trademarks or any other materials, except as permitted herein, is expressly prohibited and may be in violation of federal or other applicable laws.

© Copyright 2026 Bionano Genomics, Inc. All rights reserved.