

Bundle Component Information for 90058

First Gen Saphyr Instrument Re-Qualification Kit

Part #		Description	Qty
1	20319	Saphyr Chip, G1.2	1 ea
2	20400	DLS Biological Control A, 100µL	1 ea

Part #	Description
20400	DLS Biological Control A, 100μL
Length (cm)	1.041
Height (cm)	4.44
Volume (mL)	3.78
Weight (gram)	10
Storage Temps	+2°C +8°C
HS Code	3822.90.0000



Part #	Description
20319	Saphyr Chip, G1.2
Length (cm)	130175
Width (cm)	6.35
Height (cm)	15.875
Volume (cm³)	1312.25
Weight (gram)	99.79
Storage Temps	+2°C +8°C
HS Code	9011.90.0000

bionano™

bionano®
GENOMICS

saphyr chip®
G1.2

PN 20319
LN 18301
Exp 2024-03-02



SN VU5F-LZOL-PQO4-RN6U

Store at 2-8°C

*For Research Use Only
Not for use in diagnostic procedures*

