

Bundle Component Information for 80062

Bionano Prep SP-G2 Bone Marrow Aspirate Add-On, 12 Rxn

	Part #	Description	Qty
1	00006	Bionano Prep SP-G2 BMA Add-On (12 Reactions) Kit Insert	1 ea
2	10173	Label, Avery 5163	1 ea
3	10249	Kit Box Reagent Rack #6	1 ea
4	10307	SP Kit Box 2, 7.5" x 5.5" x 4.5"	1 ea
5	20398	DNA Stabilizer, 4 ml	1 ea
6	20452	Microcentrifuge Tubes, 2 mL, 12 ea	2 ea
7	20464	Bone Marrow Aspirate (BMA) Filters, 24 ea	1 ea
8	31002	Label,Caution Store at 15C to 30C	1 ea

Part #	Description
80062	Bionano Prep SP-G2 BMA Add-On, 12 Rxn
Length (cm)	19.675
Width (cm)	11.43
Height (cm)	13.6525
Volume (cm ³)	3070.25
Weight (gram)	172.9
Storage Temps	+15°C +25°C
HS Code	3822.90.0000

Bionano Prep SP-G2 BMA Add-On (12 Reactions)

PN 80062

LN 100000-0000000

Exp YYYY-MM-DD

Store at 15°C to 30°C

For research use only
Not for use in diagnostic procedures



QTY: 1
Box 1 of 1

Bionano Prep SP-G2 BMA Add-On (12 Reactions)

PN 80062

LN 100000-0000000

Exp YYYY-MM-DD

Store at 15°C to 30°C

For research use only
Not for use in diagnostic procedures



QTY: 1
Box 1 of 1

