

Bundle Component Information for 90215A

Stratys Core G2 Blood and Cell 48 Genome Bundle

	Part #	Description	Qty
1	80104	Stratys Core Chip, 12-pack	4 ea
2	80117-1	Bionano Prep DLS-G2 Labeling Kit, 24 Rxn, Box 1	2 ea
3	80117-2	Bionano Prep DLS-G2 Labeling Kit, 24 Rxn, Box 2	2 ea
4	80117-3	Bionano Prep DLS-G2 Labeling Kit, 24 Rxn, Box 3	2 ea
5	80118-1	Bionano Prep SP-G2 Blood & Cell Isolation Kit, 24 Rxn, Box 1	2 ea
6	80118-2	Bionano Prep SP-G2 Blood & Cell Isolation Kit, 24 Rxn, Box 2	2 ea
7	80118-3	Bionano Prep SP-G2 Blood & Cell Isolation Kit, 24 Rxn, Box 3	2 ea

Part #	Description
80104	Stratys Core Chip, 12-pack
Length (cm)	19.685
Width (cm)	11.43
Height (cm)	13.97
Volume (cm ³)	3143.24
Weight (gram)	327
Storage Temps	+2°C +8°C
HS Code	9011.90.0000

Stratys Core Chip, 12-Pack

PN 80104

SN 101010-000123401

LN 101010-0001234

Exp YYYY-MM-DD

Store at 2 to 8°C

For research use only
Not for use in diagnostic procedures



QTY: 1

Stratys Core Chip, 12-Pack

PN 80104

SN 101010-000123401

LN 101010-0001234


Exp YYYY-MM-DD

Store at 2 to 8°C

For research use only
Not for use in diagnostic procedures



QTY: 1



Part #	Description
80118-1	Bionano Prep SP-G2 Blood & Cell Isolation Kit, 24 Rxn, Box 1
Length (cm)	19.685
Width (cm)	11.43
Height (cm)	13.97
Volume (cm ³)	3143.24
Weight (gram)	478.87
Storage Temps	+15°C +25°C
HS Code	3822.90.0000



Bionano Prep SP-G2
Blood and Cell DNA
Isolation Kit (24 Reactions)

PN 80118-1
LN 100000-0000000
Exp YYYY-MM-DD



Store at 15°C to 30°C

For research use only
Not for use in diagnostic procedures

QTY: 1
Box 1 of 3

Bionano Prep SP-G2
Blood and Cell DNA
Isolation Kit (24 Reactions)

PN 80118-1
LN 100000-0000000
Exp YYYY-MM-DD



Store at 15°C to 30°C

For research use only
Not for use in diagnostic procedures

QTY: 1
Box 1 of 3

Part #	Description
80118-2	Bionano Prep SP-G2 Blood & Cell Isolation Kit, 24 Rxn, Box 1
Length (cm)	6.35
Width (cm)	5.08
Height (cm)	5.08
Volume (cm ³)	163.87
Weight (gram)	24.8
Storage Temps	+2°C +8°C
HS Code	3822.90.0000



Bionano Prep SP-G2 Blood
and Cell DNA Isolation Kit
(24 Reactions)

PN 80118-2
LN 100000-0000000
Exp YYYY-MM-DD
Store at 2°C to 8°C

For research use only
Not for use in diagnostic procedures




QTY: 1
Box 2 of 3

Bionano Prep SP-G2 Blood
and Cell DNA Isolation Kit
(24 Reactions)


PN 80118-2
LN 100000-0000000
Exp YYYY-MM-DD
Store at 2°C to 8°C

For research use only
Not for use in diagnostic procedures



QTY: 1
Box 2 of 3

Part #	Description
80118-3	Bionano Prep SP-G2 Blood & Cell Isolation Kit, 24 Rxn, Box 1
Length (cm)	6.35
Width (cm)	5.08
Height (cm)	5.08
Volume (cm ³)	163.87
Weight (gram)	24.6
Storage Temps	-20°C -25°C
HS Code	3822.90.0000



Bionano Prep SP-G2 Blood and Cell DNA Isolation Kit (24 Reactions)


PN 80118-3

LN 100000-0000000

Exp YYYY-MM-DD

Store at -25°C to -15°C

For research use only
Not for use in diagnostic procedures



QTY: 1
Box 3 of 3

Bionano Prep SP-G2 Blood and Cell DNA Isolation Kit (24 Reactions)

PN 80118-3

LN 100000-0000000

Exp YYYY-MM-DD

Store at -25°C to -15°C

For research use only
Not for use in diagnostic procedures



QTY: 1
Box 3 of 3

Part #	Description
80117-1	Bionano Prep DLS-G2 Labeling Kit, 24 Rxn, Box 1
Length (cm)	12.7
Width (cm)	6.35
Height (cm)	15.875
Volume (cm ³)	128.24
Weight (gram)	175
Storage Temps	+15°C +25°C
HS Code	3822.90.0000

Bionano Prep DLS-G2
Labeling Kit (24 Reactions)

PN 80117-1

LN XXXXXXXXX

Exp YYYY-MM-DD

Store at 15°C to 30°C

For research use only
Not for use in diagnostic procedures



QTY: 1
Box 1 of 3

Bionano Prep DLS-G2
Labeling Kit (24 Reactions)

PN 80117-1

LN XXXXXXXXX

Exp YYYY-MM-DD

Store at 15°C to 30°C

For research use only
Not for use in diagnostic procedures



QTY: 1
Box 1 of 3



Part #	Description
80117-2	Bionano Prep DLS-G2 Labeling Kit, 24 Rxn, Box 2
Length (cm)	9
Width (cm)	4.5
Height (cm)	5.1
Volume (cm ³)	206.55
Weight (gram)	28
Storage Temps	+2°C +8°C
HS Code	3822.90.0000

**Bionano Prep DLS-G2
Labeling Kit (24 Reactions)**

PN 80117-2
LN XXXXXXXXXX
Exp YYYY-MM-DD
Store at 2°C to 8°C




QTY: 1
Box 2 of 3

For research use only
Not for use in diagnostic procedures

**Bionano Prep DLS-G2
Labeling Kit (24 Reactions)**

PN 80117-2
LN XXXXXXXXXX
Exp YYYY-MM-DD
Store at 2°C to 8°C



QTY: 1
Box 2 of 3

For research use only
Not for use in diagnostic procedures



Part #	Description
80117-3	Bionano Prep DLS-G2 Labeling Kit, 24 Rxn, Box 3
Length (cm)	9
Width (cm)	4.5
Height (cm)	5.1
Volume (cm ³)	206.55
Weight (gram)	31
Storage Temps	-20°C -25°C
HS Code	3822.90.0000

Bionano Prep DLS-G2
Labeling Kit (24Reactions)

PN 80117-3

LN XXXXXXXXXX

Exp YYYY-MM-DD

Store at -25°C to -15°C

For research use only

Not for use in diagnostic procedures



QTY: 1
Box 3 of 3

Bionano Prep DLS-G2
Labeling Kit (24Reactions)

PN 80117-3

LN XXXXXXXXXX

Exp YYYY-MM-DD

Store at -25°C to -15°C

For research use only

Not for use in diagnostic procedures



QTY: 1
Box 3 of 3

